

## Controlled Heating Mantles (E Series)

**Brand:** Medline Scientific

Energy-regulated, all-purpose heating mantles.

Their time proportional double pole systems allow stable control of the heating element, reducing the risk of overheating liquids. The dial controller allows accurate temperature settings, Retort rod holders are integrated into the base of the units to accommodate our range of accessories.



Specifications:

Product Code(s)	E-1005N	E-101N	E-102N	E-103N	E-104N	E-105N	E-106N	E-107N	E-108N	E-109N
Case Construction	Steel, with chemical resistant coating									
Thermal Insulation	Moulded silica powder									
Heating Element	Nichrome wire									
Mantle	Woven fibre glass mat									
Element Temperature	Up to 450°C									
Power Supply	220V AC									
Volume (ml)	50	100	250	500	1000	2000	3000	5000	10,000	20,000
Watts	60 W	80 W	150 W	260 W	350 W	450 W	550 W	700 W	1200 W	1800 W



# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

