


**LE410****3-in-1****pH****Order number** 51340348

former order number -

**Specifications**

pH range 0 - 14

zero point (pH) 7.00 ± 0.25

temperature range 0 - 80 °C

temperature probe NTC30K

shaft material Glass

membrane glass HA

membrane resistance &lt; 600 MΩ

reference system Ag / AgCl

reference diaphragm Ceramic

reference electrolyte 3mol/L KCl

AgCl saturated

bridge electrolyte -

electrode head LE cap

cable length 1.0 m

connector BNC / Cinch

cell constant (cm<sup>-1</sup>)

cell type -

indicating metal -

**Storage**

3 mol/L KCl

**Special features**

Fast measurement in aqueous solution

Integrated temperature probe

Refillable reference electrolyte

Scale 1:1



# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

