

**Miele**

# 120 years of perfection.

**Miele Professional. Immer Besser.**

SERIES 120



Celebrate with us and profit  
from beneficial discounts on  
Miele lab washers\*

\* The campaign covers models PG 8581 and PG 8591  
and is due to run for 120 days. Campaign period: March 1, 2019 until June 29, 2019.  
For more information: 0800 22 44 644, [www.miele.de](http://www.miele.de)



## Tradition, quality and durability - over 120 years

Miele Professional develops commercial solutions to reduce your everyday workload. Our solutions are designed to meet the highest of demands, and have done so now for the past 120 years. Miele stands for uncompromising quality and durability. One thing is always at the centre of Miele's attention, at all times and no matter what we are doing: The reliable reprocessing of laboratory glassware in readiness for analytical experiments. Leaving you to concentrate on your core business.

The Miele 'Series 120' anniversary model can be combined with the innovative EasyLoad basket system and hence sets new standards in terms of convenience and efficiency. The new injector system with its innovative plastic supports enable the highly convenient positioning even of large-volume laboratory glasses whilst utilising the space in the chamber to the full.

\*Campaign period: March 1, 2019 – June 30, 2019. This campaign offer applies only to commercial/self-employed clients purchasing during the campaign period. This campaign only applies to machines purchased in Germany.



Miele lab washers



EasyLoad system

## Sample combinations: \*

- 'Series 120' campaign model
- 1 x A 101 upper basket, height-adjustable, with integrated spray arm
- 1 x A 150 lower basket, height-adjustable, for connection of injector modules
- 2 x A 300/2 EasyLoad injector modules

\* The accessories listed below do not constitute part of the 'Series 120' promotional model.

- 2 x A 860 EasyLoad frames
- 8 EasyLoad injector nozzles, Ø 6 x 130 mm
- 8 EasyLoad injector nozzles, Ø 6 x 210 mm
- A 802 nozzle, for optimum rinsing of dispenser compartment
- 2 x AK 12/1 inserts, for various utensils

Note: 'Series 120' lab washers are not suitable for reprocessing pipettes, laboratory glassware contaminated with oil or for applications requiring a temperature in excess of 70°C, e.g. to remove agar.



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

