

Connection and Commissioning for Miele Glass Washers

To connect to existing services as the terms and conditions below and carry out commissioning to the following models G7883, G7893, G7883CD, PG8535

To connect to existing services as the terms and conditions below and carry out commissioning to the following model G7804

To connect to existing services as the terms and conditions below install a plinth and carry out commissioning to all the models above

Terms and conditions

All connections to the mains electricity supply either 230 volt or 415 volt will be made with IP44 "commando type socket" (see note 1)

All hot and cold water connections are to have control valves with 3/4 BSP male fittings and be within one meter.

The deionised water connection is adjacent to the machine within 1.5 meters.

The new machine is unboxed and placed in the vicinity of the final position.

For 900mm machines to be installed on a plinth assistance would be required this is a two person operation.

Notes

1, If replacing an existing appliance that has been hard wired it is acceptable (as long as the old lead is of good order) for it to be cut and an IP44 socket to be fitted in line.

2, For a G7804 just a 13 amp socket is required.

3, If assistance cannot be offered to install a 900mm plinth an extra charge will be levied for an extra Miele technician to attend.