


CAT. NO.....	169-0045		167-0045	
PORE SIZE (µm).....	0.45		0.45	
UPPER CAP. (mL).....	500		1000	
RECEIVER CAP. (mL)...	500		1000	
MEMBRANE COLOR.....	WHITE		WHITE	
MEMBRANE DIA (mm)...	90		90	
	INCHES	MM	INCHES	MM
A (UPPER I.D.).....	3.84	97.5	4.48	113.8
B HEIGHT.....	10.7	272	12.8	325
C DIA (BASE).....	3.86	98.0	4.48	113.8
D OVERALL.....	4.0	102	4.2	106
NOM. WT. (gms).....	170		250	
MATERIALS				
UPPER CONTAINER	PS			
COVER	PS			
COLLAR	BLUE PS			
RECEIVER	PS			
STORAGE CLOSURE	WHITE HDPE			
ADAPTOR	WHITE LDPE			
GASKET	FOAMED LDPE			
MEMBRANE	PES			
VENT PLUG	CA			

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	MC ²	APP'D	PVC	PSS 00-016X-0045-4	REV. 1205
CK'D	SH		12/22/05		
DESCRIPTION: PES FILTER UNITS					



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

