


CAT. NO.....:	565-0020		568-0020		566-0020	
UPPER CAP. (mL)....:	150		250		500	
RECEIVER CAP. (mL)..:	150		250		500	
PORE SIZE (µm).....:	0.2		0.2		0.2	
MEMBRANE COLOR.....:	WHITE		WHITE		WHITE	
MEMBRANE DIA (mm)...:	50		50		75	
	INCHES	MM	INCHES	MM	INCHES	MM
A (UPPER I.D.).....:	3.15	80.0	3.15	80.0	3.80	96.5
B HEIGHT.....:	6.15	156.2	7.9	200	9.8	249
C DIA (BASE).....:	3.50	88.9	3.50	88.9	3.86	98.0
D OVERALL.....:	3.7	94	3.7	94	4.0	102
NOM. WT. (gms).....:	78		95		155	

MATERIALS	
UPPER CONTAINER	PS
COVER	PS
COLLAR	BLUE PS
RECEIVER	PS
STORAGE CLOSURE	WHITE HDPE
ADAPTOR	WHITE LDPE
GASKET	FOAMED LDPE
MEMBRANE	αPES
VENT PLUG	CA

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	PB	APP'D	JMC	PSS 00-056X-0020-1	REV.	1107
CK'D	JMC		11/29/07			
DESCRIPTION: αPES FILTER UNITS						



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

