


CATALOG NUMBER	597-3320		597-4520	
PORE SIZE (um)	0.2		0.2	
UPPER CAP. (mL)	1000		1000	
BOTTLE NECK SIZE	33mm		45mm	
MEMBRANE COLOR	WHITE		WHITE	
MEMBRANE DIA. (mm)	90		90	
	INCH	MM	INCH	MM
A (UPPER I.D.).....:	4.49	114.0	4.49	114.0
B (HEIGHT).....:	6.43	163.3	6.52	165.6
C (DIA.).....:	4.80	121.9	4.80	121.9
D OVERALL.....:	4.22	107.2	4.48	113.8
NOM. WT. (gms).....:	118		120	
	MATERIALS			
UPPER CONTAINER	PS		PS	
COVER	PS		PS	
COLLAR	BLUE PS		BLUE PS	
ADAPTOR	WHITE LDPE		WHITE LDPE	
GASKET	FOAMED LDPE		FOAMED LDPE	
MEMBRANE	PES		PES	
VENT PLUG	CA		CA	

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	PB	APP'D	GRP	PSS 00-0597-XX20	REV. 0613
CK'D	KOW	6/26/13			
DESCRIPTION:				BOTTLE TOP PES FILTER	



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

