

# Chest freezer LT C300

- Low cost of ownership
- Limited maintenance
- On-site service
- Easy removable

Specifications	
<b>Technical</b>	
Gross volume:	296 litres
Temp. range:	-10°C to -45°C
Refrigerant:	R404a
Baskets:	1
Power consumption 24/h:	5,3 kWh
Voltage:	220V-240V 50/60Hz
Compressor type:	Single stage
<b>Features</b>	
- Audible and visual alarm	
- High/low alarm	
- Battery backup (optional)	
- Prepared for CO2 backup	
- Castors	
- Handle lock	
<b>Outer dimensions</b>	
Height:	885 mm
Width:	1260 mm
Depth:	695 mm
Gross weight:	82 kg
Net weight:	76 kg
<b>Inner dimensions</b>	
Height:	635/385 mm
Width:	890/210 mm
Depth:	440 mm
Product no.	N111003



Storage capacity				
Box height	No. of boxes	No. of racks	Total no. of boxes	Product no.
50 mm	6	3	18	N311200
	11	18	198	N311220
			<b>= 216</b>	
75 mm	4	3	12	N311300
	8	18	144	N311320
			<b>= 156</b>	
100 mm	3	3	9	N311400
	5	18	90	N311420
			<b>= 99</b>	



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

