

Chest freezer LT C400

- Low cost of ownership
- Limited maintenance
- On-site service
- Easy removable

Specifications	
Technical	
Gross volume:	383 litres
Temp. range:	-10°C to -45°C
Refrigerant:	R404a
Baskets:	1
Power consumption 24/h:	5,8 kWh
Voltage:	220V-240V 50/60Hz
Compressor type:	Single stage
Features	
- Audible and visual alarm	
- High/low alarm	
- Battery backup (optional)	
- Prepared for CO2 backup	
- Castors	
- Handle lock	
Outer dimensions	
Height:	885 mm
Width:	1560 mm
Depth:	695 mm
Gross weight:	95 kg
Net weight:	88 kg
Inner dimensions	
Height:	635/385 mm
Width:	1190/210 mm
Depth:	440 mm
Product no.	N111004



Storage capacity				
Box height	No. of boxes	No. of racks	Total no. of boxes	Product no.
50 mm	6	3	18	N311200
	11	24	264	N311220
			= 282	
75 mm	4	3	12	N311300
	8	24	192	N311320
			= 204	
100 mm	3	3	9	N311400
	5	24	120	N311420
			= 129	



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

