

Chest freezer LT C500

- Low cost of ownership
- Limited maintenance
- On-site service
- Easy removable

Specifications	
Technical	
Gross volume:	495 litres
Temp. range:	-10°C to -45°C
Refrigerant:	R404a
Baskets:	3
Power consumption 24/h:	6,5 kWh
Voltage:	220V-240V 50/60Hz
Compressor type:	Single stage
Features	
- Audible and visual alarm	
- High/low alarm	
- Battery backup (optional)	
- Prepared for CO2 backup	
- Legs	
- Handle lock	
Outer dimensions	
Height:	890 mm
Width:	1660 mm
Depth:	758 mm
Gross weight:	122 kg
Net weight:	107 kg
Inner dimensions	
Height:	640/390 mm
Width:	1290/210 mm
Depth:	500 mm
Product no.	N111005



Storage capacity				
Box height	No. of boxes	No. of racks	Total no. of boxes	Product no.
50 mm	6	3	18	N311200
	11	27	297	N311220
			= 315	
75 mm	4	3	12	N311300
	8	27	216	N311320
			= 228	
100 mm	3	3	9	N311400
	5	27	135	N311420
			= 144	



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

