

Upright freezer LT U250

- Low cost of ownership
- Limited maintenance
- On-site service
- Easy removable

Specifications	
Technical	
Gross volume:	253 litres
Temp. range:	-10°C to -45°C
Refrigerant:	Nature R
Inner doors:	-
Adjustable shelves:	5
Power consumption 24/h:	5,5 kWh
Voltage:	220V-240V 50/60Hz
Compressor type:	Single stage
Features	
- Temperature logger with USB port	
- Door alarm	
- High/low alarm	
- 48h battery backup (optional)	
- Castors	
Outer dimensions	
Height:	2055 mm
Width:	600 mm
Depth:	595 mm
Gross weight:	125 kg
Net weight:	112 kg
Inner dimensions	
Height:	1420 mm
Width:	420 mm
Depth:	370 mm
Product no.	N121002



Storage capacity				
Box height	No. of boxes	No. of racks	Total no. of boxes	Product no.
50 mm	4 x 2 = 8	3	24	N321200
	5 x 2 = 10	12	120	N321210
			= 144	
75 mm	2 x 2 = 4	3	12	N321300
	3 x 2 = 6	12	72	N321310
			= 84	
100 mm	2 x 2 = 4	3	12	N321400
	2 x 2 = 4	12	48	N321400
			= 60	



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

