

Mini Hotplates & Stirrers



OHAUS Mini hotplates and stirrers are rugged, compact units that can heat and stir up to 1000 ml of water. These are ideal for educational and other lab environments that require consistent heating and stirring. Powerful heater reaches maximum temperature in minutes. Models include Hotplates, Stirrers, Hotplate-Stirrers, Fixed Temperature Hotplates and Auto Stirrers. All models feature a 10.2 × 10.2 cm, easy-to-clean ceramic top plate and cool-touch housing.

- Chemical-Resistant, Easy-to-Clean Ceramic Top Plate
- Constructed of Cool-Touch, Chemical-Resistant Housing
- Equipped with a Built-In Support Rod Holder

Hotplates & Stirrers

Mini Hotplates & Stirrers

- New ergonomic design
- Hotplates & Hotplate-Stirrers boil 300 mL of water in 18 minutes
- Ideal for educational labs
- Built in support rod holder

The OHAUS Basic Mini Hotplates, Stirrers, and Hotplate-Stirrers are rugged, compact units that heat and stir up to 1000 mL of liquid. Ceramic tops feature a chemical resistant white reflective top plate that is easy to clean. Bi-metallic thermostat offers reliable temperature control. Powerful heater reaches maximum temperature in only minutes. Powerful motor and magnet deliver reliable and consistent stirring. Compact design saves bench space. Built-in support rod holder with locking knob accepts optional Support Rod and Clamp Kit.

Basic Mini Fixed Temperature Hotplate features an illuminated rocker switch to activate the preset fixed temperature of 500°C.

Basic Mini Auto-Stirrer is automatically activated by the minimum weight of a beaker or flask and will stop stirring when mixing vessel is removed.

Operating Features:

Adjustment Knobs: Basic speed and temperature control knobs with dial markings from 1 to 10.

Operating Conditions:

Units can be run in conditions from 5 to 40°C, 20% to 85% relative humidity, non-condensing.

Applications:

Academia and basic chemistry.

Ordering Information:

Units include a 3-wire cord and plug. Stirrers and Hotplate-Stirrers are supplied with a 3.8 cm PTFE coated stir bar.

Basic Mini Hotplates

Description	Model	Item Number
Basic Mini Hotplate	HSMNHP4CAL	30392012
Basic Mini Fixed Temperature Hotplate	HSMNHP4CFT	30392033

Basic Mini Stirrers

Description	Model	Item Number
Basic Mini Stirrer	HSMNST4CAL	30392019
Basic Mini Auto-Stirrer	HSMNAS4CAL	30392040

Basic Mini Hotplate-Stirrers

Description	Model	Item Number
Basic Mini Hotplate-Stirrer	HSMNHS4CAL	30392026



Specifications	
Temperature Range	to 500°C*
Speed Range	100 to 1200 rpm
Maximum Capacity (H ₂ O)	1000 mL
Top Plate Dimensions (Dia.)	10.2 × 10.2 cm
Overall Dimensions (L × W × H)	20.1 × 15.2 × 12.7 cm
Ship Weight	1.8 kg

* Fixed Temperature Hotplate has a fixed temperature of 500°C.





Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

