



Model:
SFX-34

Typology:
EDUCATIONAL STEREOMICROSCOPE

Description:
Entry- level stereomicroscope for routine and educational fields.

Head	Binocular, 45° inclined. Adjustable dioptric compensation on both eyepiece tubes. Interpupillary distance adjustable between 53 and 75 mm.
Eyepieces	Wide Field, WF 10x/20.
Objectives	1x – 3x selectable. Working distance: 57 mm.
Stand	Fixed with focus mechanism, with handle.
Focusing	By means of a rack and pinion controlled by a pair of knobs placed on both sides of the microscope.
Stage	Fitted with a semi-transparent plastic disc for transmitted light with specimen clamps.
Illumination	Double adjustable illuminator for incident and transmitted light fitted with two high-efficiency 1 W LED light sources. Touch brightness control. Input Voltage: 110/240 Vac, 50/60 Hz. Maximum power consumption: 3 W. Internal NiMH rechargeable battery. Up to 8 hours of use from a single charge, LED at medium intensity. Recharging time: approx. 8 hours.
Dimensions	HEIGHT: 288 mm WIDTH: 183 mm DEPTH: 205 mm WEIGHT: 1.8 kg
Accessories included	Manual and dust cover.



Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

