



Model:  
**SFX-52**

Typology:  
**EDUCATIONAL STEREO MICROSCOPE**

**Description:**  
*Entry- level stereomicroscope for routine and educational fields.*

<b>Head</b>	Binocular, 45° inclined, 360° rotating. Adjustable dioptric compensation on both eyepiece tubes. Interpupillary distance adjustable between 53 and 75 mm.
<b>Eyepieces</b>	Wide Field, WF 10x/20.
<b>Objectives</b>	1x – 3x selectable. Working distance: 76 mm.
<b>Stand</b>	Fixed with focus mechanism, with handle.
<b>Focusing</b>	By means of a rack and pinion controlled by a pair of knobs placed on both sides of the microscope.
<b>Stage</b>	Fitted with a semi-transparent plastic disc for transmitted light with specimen clamps.
<b>Illumination</b>	Double adjustable illuminator for incident and transmitted light fitted with two high-efficiency 1 W LED light sources. Touch brightness control. Input Voltage: 110/240 Vac, 50/60 Hz. Maximum power consumption: 3 W. Internal NiMH rechargeable battery. Up to 8 hours of use from a single charge, LED at medium intensity. Recharging time: approx. 8 hours.
<b>Dimensions</b>	HEIGHT: 288 mm WIDTH: 183 mm DEPTH: 205 mm WEIGHT: 1.9 kg
<b>Accessories included</b>	Manual and dust cover.



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax:01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

