



Binocular stereozoom microscope,
6.7x...45x, pillar stand

| | | |
|-------------------------------------|--|--|
| Observation Method - Incident Light | Brightfield | Yes |
| | Simple polarized light | As optional |
| Head | Type | Binocular |
| | Construction material | Plastic mold / Aluminum gears |
| | Inclination | 45° |
| | 360° rotating | Yes |
| | Interpupillary distance (mm) | 51-75 |
| | Fixing screw for eyepieces | Yes |
| | Head holder internal diameter (mm) | 76 |
| | Tube inner diameter (mm) | 30 |
| Eyepieces | Field number (mm) | 23 |
| | Magnification | 10x |
| | Planar type | Yes |
| | Micrometric scale | As optional |
| | Diameter of micrometer glass (mm) | 26 |
| | High eyepoint (for glass wearers) | Yes |
| | Dioptric adjustment | Yes |
| | Rubber cup | Yes |
| | Retractable protections | Yes |
| Objective | Optical system | Greenough |
| | Anti-fungus treatment | Yes |
| | Zoom Type | Parfocal achromatic |
| | Working distance (mm) | 110 |
| | Standard magnifications | 6.7x-45x |
| | Zoom ratio | 6.72 |
| | Zoom click stops | At 0.67x, 0.8x, 1x, 1.2x, 1.5x, 2x, 2.5x, 3x, 3.5x, 4x, 4.5x |
| Focusing System | Type | Coarse |
| | Coarse total travel (mm) | 50 |
| | Adjustable tension | Yes |
| | Head holder internal diameter (mm) | 76 |
| Stand | Type | Pillar |
| | Pillar diameter (mm) | 32 |
| Stage | Black/white contrast plate for incident ill. | Yes, diameter 100 mm |
| | Specimen clips | Yes |
| Accessories Included | Dust cover | Yes |
| | User Manual | Digital version (downloadable) |
| Additional Information | Maximum sample height (mm): | 110 |
| Product Dimensions | Total height (mm) | 315 |
| | Head height (mm) | 240 |
| | Head width (mm) | 167 |
| | Head depth (mm) | 200 |
| | Base height (mm) | 30 |
| | Base width (mm) | 210 |
| | Base depth (mm) | 270 |
| Product Weight | (kg) | 4 |



Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

