



Model:

B-500POL

Typology:

RESEARCH MICROSCOPE

Description:

Laboratory microscope for research applications.

Dye-cast frame, with high stability and ergonomy, for transmitted light observation.

| | |
|--------------------------|--|
| Illumination | Light source type X-LED with white LED; light intensity control using a knob on left side of the frame. LED power 3.6W, comparable to an halogen bulb 35W. Color temperature: 6300K LED average lifetime approx.: 50.000h. The light exit can be used as a filter holder for additional filters (blue, yellow, frosted). Voltage: 110/230Vac, 50/60Hz, 0.4/0.8A; Fuse: T3.15A 250V Max power required: 7W |
| Observation Modes | Brightfield, Polarized transmitted light |
| Focusing | Coaxial coarse and fine focusing mechanism (graduated, 0.002mm) with upper stop, to prevent the contact between objective and specimen. Adjustable tension of coarse focusing knob. |
| Stage | Rotatable stage with locking system and Vernier scale with accuracy 0,1 mm. Diameter 160 mm. Specimen clamps. |
| Nosepiece | Quadruple revolving nosepiece, rotation on ball bearings. Centering system on each objective. |
| Head | Trinocular observation head, inclined 30° and rotatable 360°. Diopter adjustment on left eyepiece sleeve. Interpupillary adjustment 55-75 mm. |
| Eyepieces | Wide field eyepieces WF10X/22 with field number 22. One eyepiece has a built-in crosshair reticle for the centering of the optical path. |
| Objectives | Infinity corrected optical system IOS (Infinity Optical System). Plan-achromatic "strain-free" objectives infinity corrected, made by following objectives: -) PLAN IOS POL 4X , A.N. 0.10, W.D. 4,7 mm -) PLAN IOS POL 10X, A.N. 0.25, W.D. 4,1 mm -) PLAN IOS POL 40X, A.N. 0.65, W.D. 0,5 mm -) PLAN IOS POL 60X, A.N. 0.85, W.D. 0,1 mm All objectives are treated with an anti-fungus treatment. |
| Condenser | Abbe condenser, N.A. 0,9 with centering system. Rotatable polarizer with graduated scale (it can be excluded from optical path through rotation). |
| Bertrand Lens | Swing-out type with centering and focusing mechanism for observation in conoscopy/orthoscopy. Rotatable analyzer from 0° to 90° with graduated scale. Tint plates included: -) 1° order red (1 lambda) -) ¼ lambda -) Quartz wedge |
| Dimensions | HEIGHT: 480 mm WIDTH: 245 mm DEPTH: 380 mm WEIGHT: 9 kg |
| Accessories | Instruction manual and dust cover included. |



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

