



Model:

**POL-X**

Typology:

**LED POLARIMETER**

Description:

*Bench polarimeter with monochromatic LED light source*

Measuring range of optical rotation	$\pm 180^\circ$
Resolution	$1^\circ$
Accuracy	$0.05^\circ$
Magnification factor of the magnifying glass	4X
Light source	Monochromatic LED, 1.2 W, $\lambda_d = 590$ nm (equivalent to sodium lamp)
Length of test tubes	100mm and 200mm
Power Supply	Input 100/240Vac, 50/60 Hz Output 5Vdc 500mA
Weight	1,7kg
Overall dimensions	LxWxH = 450 x 180 x 320 mm



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

