



CORK BORERS, SPATULAS, TONGS & RACKS

SECTION 06 CORK BORERS



SPATULAS



1. Manufactured from brass and nickel plated. The sets come complete with push rod for clearing borers. Sizes numbered consecutively from 1 upwards in each set.

Catalogue Number.	Borers In Set.
3400	6
3401	12
3402	18

All sizes are 100 mm in length with a tube wall thickness of 0.45 mm



2.

2. **3403 | CORK BORER SHARPENER**
Manufactured out of aluminium with a blue handle and a hardened steel blade, for borer diameters 4 to 27 mm.

3. **3410 | CORK BORING MACHINE**
For boring holes 3 to 25 mm diameter in cork or rubber stoppers. Manufactured out of cast iron and painted blue acrylic. The borer has a 90 mm rubber base to support stoppers. The shaft has a central hole, which accommodates a push rod for clearing borers. Complete with set of 15 cork borers, sizes range from 4.76 to 27.00 mm, wall thickness 0.9 mm.

3411
Spare set of 15 cork borers and push rod. Sizes as supplied with cork boring machine 3410.

Size	Outside Dia. (mm).	Size	Outside Dia. (mm).	Size	Outside Dia. (mm).
1	5.00	7	12.50	13	20.00
2	6.25	8	13.75	14	21.25
3	7.50	9	15.00	15	22.50
4	8.75	10	16.25	16	23.75
5	10.00	11	17.50	17	25.00
6	11.25	12	18.75	18	26.25



3.



1.



5, 4.



2.



3.



6.

1. **SPATULAS 5500 SERIES.**
Manufactured in stainless steel - chattaway pattern, one end flat, the other bent.

Catalogue Number.	Blade Width.	Overall Length.
5500/1	5 mm	100 mm
5500/2	5 mm	230 mm
5500/3	8 mm	150 mm
5500/4	10 mm	200 mm

2. **SPATULAS 5501 SERIES.**
Stainless steel, micro chattaway pattern.

Catalogue Number.	Blade Width.	Overall Length.
5501/1	4 mm	100 mm
5501/2	4 mm	150 mm

3. **SPATULAS 5502 SERIES.**
Stainless steel. Spoon one end, flat the other.

Catalogue Number.	Spoon Diameter.	Overall Length.
5502/1	16 mm	150 mm
5502/2	18 mm	180 mm
5503/3	18 mm	210 mm

4. **SPATULAS 5503.**
Stainless steel. Offset spoon one end, flat the other. Spoon diameter 20 mm, overall length 190 mm.

5. **SPATULAS 5505.**
Stainless steel. Trowel pattern. Overall length 175 mm, width large end 19 mm, width narrow end 10 mm.

6. **SPATULAS 5504.**
Stainless steel. Nuffield pattern. Scoop at one end, narrow bent blade at the other. Width 10 mm, overall length 140 mm.



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

