



## SLS Select Indicating Tapes



Autoclavable tape indicates when steam sterilisation has occurred by means of a visible black strip which shows after autoclaving for 3 minutes at 134°C, or 15 minutes at 121°C.

Code	Width, mm	Length, m	Price
<b>SLS1610</b>	19	50	£8.40
<b>SLS1612</b>	25	50	£11.10



## 3M™ Autoclave Indicating Tapes



Specially prepared to react only under autoclave pressure conditions. A black indicator strip appears only when this tape has been steam sterilised for 3 minutes at 134°C or 15 minutes at 121°C. Rolls are individually wrapped.

Code	Alt Ref	Width, mm	Length, m	Price
<b>AUT1610</b>	1322-18	18	55	£9.68
<b>AUT1612</b>	1322-24	25	55	£15.50



# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

