

## SLS SELECT Ice Bucket



Its 25mm thick insulated double wall and insert cover are made of smooth PVC inside and outside with a foam insulation sandwich between. The container has a tapered shape.

Code	Description	Int. Dims, w x d, mm	Ext. Dims, w x d, mm
ICE3200	Ice bucket	190 (top), 110 (bottom) x 150	280 x 210



## Azlon Ice Buckets



- Black, green and blue skinned polyurethane foam ice buckets
- Suitable for use with water-ice and salt solutions
- Ideal for separating work areas by colour
- Supplied with lids
- Integral handles allow easy lifting and carrying
- Octagonal shape for easy stacking and stability
- Capacity 4.5L
- Dims, h x w, mm 185 x 330
- Weight 1.2kg

**Warning! Polyurethane is not suitable for use with liquid nitrogen or solvents**

Code	Alt Ref	Description
ICE1002	IBB001P	Black
ICE1004	IBB002P	Pale blue
ICE1006	IBB003P	Green



## Ice Pack

Blue.

Code	Description
DAI1108	Ice pack



## Cool Containers



Rigid single mould polyurethane provides durability for chilling your temperature sensitive samples. Excellent insulation properties for ultra cold-use with ice, ice-salt slurries, and more. Minimise evaporation loss with snug fit lid. Keyed lid and base fit together to conveniently store lid under base. Use handy pour spout with melted ice, water, and solutions. Ergonomic handles safely transport bucket or pan from work area to sink. Tapered sides make stacking and storage compact. Suitable for use down to -196°C.

Code	Alt Ref	Style	Size, L	Colour
ICE1100	HS28714K	Bucket	4	Black Onyx
ICE1102	HS28721G	Pan	1	Green Emerald
ICE1104	HS28724B	Pan	4	Arctic Blue
ICE1106	HS28729R	Pan	9	Ruby Red





**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

**The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed, please provide these to our sales team so that details can be confirmed.**

