

SciQuip Ovens

SciQuip Ovens S-Series HT-Series



- Capacity 30ltr – 230ltr
- Temp up to 300°C (HT-Series)
- Fully programmable
- Fan-assisted excellent temp uniformity
- Audible & visual alarms
- Password protection
- Large LCD screen
- Separate overheat safety protection
- Viewing window

SciQuip fan assisted drying ovens have been designed to suit all drying applications including aging, bonding, curing, drying glassware, baking, sterilizing, heat treating, annealing, stress relieving and stability testing. The standard S-series provide temperatures up to 200°C, while the high temperature HT-series go to 300°C.

The SciQuip ovens can store up to 9 programs of 18 steps, which can be set for multiple ramp/soak cycles or for running one temperature and switching automatically to another for a given time. Each program can be made to repeat from 1 to 99 times automatically.

All S-Series ovens stand on feet, with the exception of the 210S, which is on lockable castor wheels. All HT-Series come with lockable castor wheels.



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

