

## ***SciQuip Microplate Shaker***



## SciQuip Microplate Mixer

- Can hold either one or two microtiter plates depending on platform
- Universal platform mount will accept additional accessories for shaking sample tubes
- Adapter detection
- Continuous or touch operation
- Stepless speed regulation up to a maximum of 3000rpm with speed feedback
- The small (4.5mm diameter) orbiting motion creates true vortex uniform mixing in each well of the microplate
- Brushless motor, quiet and maintenance free

### Ordering Information

Part No.	Description
SP2260-MP	SciQuip Microplate Mixer with single plate clamp

### Accessories

Part No.	Description
SP2260-MP-2	Double microplate clamp
SP2260-PS1.1	Universal round top plate for tube adaptors
SP2260-VT1.3.1	Tube adapter, for 48 holes test tubes, Ø6mm (0.5ml)
SP2260-VT1.3.2	Tube adapter, for 18 holes test tubes, Ø10mm (1.5-2.0ml)
SP2260-VT1.3.3	Tube adapter, for 12 holes test tubes, Ø12mm (2-7ml)
SP2260-VT1.3.4	Tube adapter, for 8 holes test tubes, Ø16mm (4-15ml)
SP2260-VT1.3.5	Tube adapter, for 8 holes test tubes, Ø20mm (20ml)

### Technical Specifications

Shaking movement	Orbital
Orbital diameter	4.5 mm
Max. shaking weight (with attachments)	0.5 kg
Speed range	0 – 3000 rpm
Speed display	Scale
Dimensions [W×D×H]	150 × 260 × 80 mm
Weight	3 kg
Power	85 - 265V 50/60Hz – 20W





**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

