



P11/ROD/DJ/BNC Specification

PH Range:	0-14
E O Point:	-25mV to +25mV
Temp Range °C:	0-80
Response Time:	<30 secs 95% FSD
Ref. Type:	AgCl, Double Junction Liquid filled
Junction Type:	Ceramic
Shaft Length:	120mm ± 2mm
Shaft Diameter:	-11.8 – 12.2mm
Cable Length:	1 metre
Connector:	BNC (others available upon request)
Combination:	Yes
Body Type:	Glass
Application:	High Accuracy measurements in Liquids and Non Aqueous samples
Refill Solution:	551-64
Storage Solution:	791-65
Recommended Cleaning Solution:	791-15
Shelf Life:	12 months



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

