



## **P11 Protected Bulb Specification**

PH Range:	0-14
E O Point;	0 to $\pm$ 25mV
Temp Range °C:	0-100
Response Time:	<30 secs 95% FSD
Ref. Type:	Single Junction Ag/AgCl
Junction Type:	Porous Teflon
Shaft Length:	120mm
Stem Dia mm:	12mm
Cable Length:	1 metre
Connector:	BNC
Combination:	Yes
Body Type:	Glass
Application:	pH measurements where a very robust device is required
Refill Solution:	551-63 - PG
Storage Solution:	791-65
Recommended Cleaning Solution:	791-15
Shelf Life:	12 months



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

