



## **P16 Specification**

PH Range:	0-14
E O Point:	-25mV to +25mV
Temp Range °C:	0-50
Response Time:	<30 secs 95% FSD
Ref. Type:	HgCl liquid filled Single Junction
Junction Type:	Frit Ceramic
Lower Shaft Length:	90mm ± 2mm
Lower Shaft Dia mm:	6mm
Cable Length:	1 Metre
Connector:	BNC (others available upon request)
Combination:	Yes
Body Type:	Glass
Application:	Tris Buffers
Refill Solution:	551-64
Storage Solution:	791-65
Recommended Cleaning Solution:	791-15
Shelf Life:	12 months



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

