



PI10/LAB Specification

pH Range:	0-14
E O Point:	-25mV to +25mV
Temp Range:	0-130°C
Ref Type:	AgCl, Double Junction
Junction Type:	Porous Teflon
Shaft Length:	120mm \pm 2mm
Shaft Dia:	11.8-12.2mm
Combination:	Yes
Cable Length:	1 Metre
Connector:	BNC (others available upon request)
Body Type:	Glass
Application:	Aggressive Samples, High pH, resists poisoning and fouling.
Storage Solution:	791-65
Recommended Cleaning Solution:	791-15
Shelf Life:	12 months



Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

