

## SciSpin Micro

- Compact design
- Maximum 15,000rpm/15,100g
- Capacity 12 x 1.5/2.2ml
- User friendly large LCD display
- Unique design minimises chamber temp
- Brushless motor
- Over speed & over temp detector
- Run-end alarm, pulse button, automatic lid open

### Ordering Information

#### SciSpin Microfuge

Part No.	Description
SQ-6000	SciSpin Micro including 12-place rotor with snap on lid

#### Accessories

Part No.	Description
SQ-6010	Rotor for 12 x 1.5/2.0ml with snap on lid
SQ-6020	0.5ml adapter for use with rotor
SQ-6025	0.2ml adapter for use with rotor

#### SciSpin Micro

Capacity	12 x 1.5/2.2ml
Max Speed	15,000rpm (increment 100rpm)
Min speed	500rpm
Max RCF	15,100g (increment 100g)
Timer	30 sec, 99 min or continuous
Dimensions (WxDxH)	255 x 245 x 140mm
Weight	6kg
Noise level	56dB
Power	220V 50/60Hz





# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

