

CRYOLAB 25 inc accessories

For dispensing into small
dewars and flasks.

The following is included:

- Cryolab 25.
(25 litres of LN₂)
- Withdrawal Device inc
Swan neck
& Phase separator.
- Roller base.

When ordering
request **Part number 9901130**



1. Easy handling of the 25 litre Dewar weighing 26 kg
2. Below PED, so no yearly safety checks required- saving money.
3. Aluminium Dewar so a lighter construction than stainless steel vessel and easier to move.
4. Can be filled by any gas company.
5. Dispense up to 2 litres/minute with limited head pressure.
6. Dispenser fitted with phase separator.

CRYOLAB 25 inc accessories

For refrigerator top up.

The following is included:

- Cryolab 25.
(25 litres of LN2)
- Withdrawal Device inc
a 1.2 Meter hose
& Phase separator.
- Roller base.
Ease of movement.

When ordering
request **Part number 9901131**



1. Easy handling of the 25 litre Dewar weighing 26kg
2. Below PED, so no yearly safety checks required- save on money.
3. Aluminium Dewar so a lighter construction than stainless steel vessel and easier to move.
4. Can be filled by any gas company.
5. Dispense up to 2 litres/minute with limited head pressure.
6. Dispenser fitted with hose and phase separator.

CRYOLAB 25 inc accessories

For refrigerator top up.

The following is included:

- Cryolab 25.
(25 litres of LN₂)
- Withdrawal Device inc
a 1.2 Meter hose
& Phase separator.
- Transport trolley.
Ease of movement.

When ordering
request **Part number 9901132**



The pouring Trolley is
1033mm high x 472mm wide.

1. Easy handling of the 25 litre Dewar weighing 26kg
2. Below PED, so no yearly safety checks required- save on money.
3. Aluminium Dewar so a lighter construction than stainless steel vessel and easier to move.
4. Can be filled by any gas company.
5. Dispense up to 2 litres/minute with limited head pressure.
6. Dispenser fitted with hose and phase separator.

CRYOLAB 25 inc accessories

For refrigerator top up.

The following is included:

- Cryolab 25.
(25 litres of LN2)
- Withdrawal Device
including Swan neck
& Phase separator.
- Transport trolley.
Ease of movement.

When ordering
request **Part number 9901133**



The pouring Trolley is
1033mm high x 472mm wide.

1. Easy handling of the 25 litre Dewar weighing 26 kg
2. Below PED, so no yearly safety checks required- save on money.
3. Aluminium Dewar so a lighter construction than stainless steel vessel and easier to move.
4. Can be filled by any gas company.
5. Dispense up to 2 litres/minute with limited head pressure.
6. Dispenser fitted with phase separator.

CRYOLAB 35 inc accessories

For dispensing into small
dewars and flasks.

The following is included:

- Cryolab 35.
(35 litres of LN₂)
- Withdrawal Device
including Swan neck
& Phase separator.
- Roller base

When ordering
request **Part number 9901134**

1. Easy to manoeuvre trolley.
2. Below PED, so no yearly safety checks required- save on money.
3. Aluminium Dewar so a lighter construction than stainless steel vessel and easier to move.
4. Can be filled by any gas company.
5. Dispense up to 2 litres/minute with limited head pressure.
6. Dispenser fitted with phase separator.

The following part numbers are included in this system:
9901080 - 9701057A - 9902082



CRYOLAB 35 inc accessories

For refrigerator top up.

The following is included:

- Cryolab 35.
(35 litres of LN₂)
- Withdrawal Device
including a
1.2 Meter hose
& Phase separator.
- Roller base

When ordering
request **Part number 9901135**

1. Easy to manoeuvre trolley.
2. Below PED, so no yearly safety checks required- save on money.
3. Aluminium Dewar so a lighter construction than stainless steel vessel and easier to move.
4. Can be filled by any gas company.
5. Dispense up to 2 litres/minute with limited head pressure.
6. Dispenser fitted with hose and phase separator.

The following part numbers are included in this system:
9901080 - 9701152A - 9902082



CRYOLAB 50 inc accessories

For dispensing into small
dewars and flasks.

The following is included:

- Cryolab 50.
(50 litres of LN2)
- Withdrawal Device
including Swan neck
& Phase separator.
- Roller base

When ordering
request **Part number 9901136**

1. Easy to manoeuvre trolley.
2. Below PED, so no yearly safety checks required- save on money.
3. Aluminium Dewar so a lighter construction than stainless steel vessel and easier to move.
4. 50% the cost of a stainless steel vessel and will keep LN2 longer.
5. Can be filled by any gas company.
6. Dispense up to 2 litres/minute with limited head pressure.
7. Dispenser fitted with phase separator.

The following part numbers are included in this system:
9901090-9701057-9902082



CRYOLAB 50 inc accessories

For refrigerator top up.

The following is included:

- Cryolab 50.
(50 litres of LN2)
- Withdrawal Device
including a
1.2 Meter hose
& Phase separator.
- Roller base

When ordering
request **Part number 9901137**

1. Easy to manoeuvre trolley.
2. Below PED, so no yearly safety checks required- save on money.
3. Aluminium Dewar so a lighter construction than stainless steel vessel and easier to move.
4. 50% the cost of a stainless steel vessel and will keep LN2 longer.
5. Can be filled by any gas company.
6. Dispense up to 2 litres/minute with limited head pressure.
7. Dispenser fitted with hose and phase separator.

The following part numbers are included in this system:
9901090-9701152-9902082





Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

