

Thermo Scientific™ Sterilin™ Serological Pipettes

Thermo
SCIENTIFIC

- Gamma irradiated to ensure sterility
- Certified non-pyrogenic with an endotoxin level <0.01EU/mL
- Batch-specific certificate supplied in each carton containing lot number for traceability
- Accuracy of $\pm 1\%$ for nominal volume
- Clear black printed graduations for maximum clarity
- Descending and ascending graduations to aid pipetting, with negative graduations for extra capacity
- Colour-coded for ease of identification
- Validated non-haemolytic



Code	Alt Ref	Description	Capacity, mL	Inner Pack	Pack	Price
PIP4001	40101K	Indiv. wrapped	1	1	1000	£288.00
PIP4003	40301K	Bulk packed	1	50	1000	£273.00
PIP4005	41301K	Open ended, bulk wrapped	1	50	1000	£281.00
PIP4007	40102K	Indiv. wrap	2	1	1000	£335.00
PIP4010	40105	Indiv. wrap	5	1	200	£107.00
PIP4012	40305	Bulk packed	5	25	200	£107.00
PIP4015	47105	Indiv. wrap, suction adapter	5	1	200	£120.00
PIP4017	47105N	Narrow orifice, suction adapter	5	50	200	£127.00
PIP4020	47110	Plastic wrap, suction adapter	10	50	500	£301.00
PIP4024	47310	Suction adapter, bulk packed	10	25	500	£299.00
PIP4029	47110N	Narrow orifice, suction adapter	10	50	500	£301.00
PIP4030	40125	Plastic wrap, suction adapter	25	25	200	£269.00
PIP4032	47150	Indiv. wrap, suction adapter	50	1	50	£184.00



Thermo Scientific™ Nunc™ Serological Pipettes

Thermo
SCIENTIFIC

- Made from 100% virgin polystyrene for maximum clarity
- Bold clear graduations with easy-to-read scale markings
- Non-dripping tip for pipette precision
- Certified non-pyrogenic
- Non-cytotoxic
- Plugged
- Sterile, gamma irradiated at a validated dose to ensure SAL of 10^{-6}
- Colour-coded packaging for ease in sorting and selecting the correct size

Code	Capacity, mL	Inner Pack	Pack	Price
170366	5	1	200	£51.00
170373	5	25	500	£125.00
170367	10	1	200	£64.00
170374	10	25	500	£150.00
170368	25	1	200	£114.00
170375	25	25	200	£114.00



Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

