

## SLS Select Serological Pipettes



- Manufactured in compliance with the current versions of the ISO 13485 and FDA CFR 21 Part 820, cGMP
- Made from USP Class VI polystyrene resin
- Crisp, dark ascending and descending as well as negative graduations
- Colour-coded rings for easy identification
- Plugged
- Non-pyrogenic
- RNase-/DNase-free
- Sterile




Code	Volume, mL	Packaging	Inner Pack	Pack	Price
SLS4054	1	Individually wrapped	100	1000	£143.00
SLS4056	1	Bulk	25	1000	£118.00
SLS4058	2	Individually wrapped	100	1000	£157.00
SLS4060	2	Bulk	25	1000	£118.00
SLS4062	5	Bulk	25	500	£112.00
SLS4064	5	Individually wrapped	50	200	£43.00
SLS4066	10	Individually wrapped	50	200	£45.00
SLS4068	10	Bulk	25	500	£137.00
SLS4070	25	Bulk	20	200	£146.00
SLS4072	25	Individually wrapped	50	200	£103.00
SLS4074	50	Individually wrapped	20	100	£126.00

## Thermo Scientific™ Sterilin™ Shortie Serological Pipettes



- Gamma irradiated to ensure sterility
- Certified non-pyrogenic with an endotoxin level <0.01EU/mL
- Batch-specific certificate supplied in each carton containing lot number for complete traceability
- Accuracy of  $\pm 1\%$  for nominal volume
- Clear black printed graduations for maximum clarity
- Descending and ascending graduations to aid pipetting, with negative graduations for extra capacity
- Colour-coded for ease of identification
- Validated non-haemolytic
- Supplied with suction adapter, unless otherwise stated

Supplied in packs of 200



Code	Alt Ref	Capacity, mL	Negative Graduations, mL	Price
PIP4016	42505	5	-3.0	£245.00
PIP4028	42510	10	-10.0	£363.00



# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

