


# Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



## Bibbypet Bottle top Dispensers BD10, BD25

- Economical dispensing of non-aggressive liquids
- Low-cost replaceable piston/cylinder and valves
- Easy cleaning and maintenance
- Simple to operate
- Automatic refill



The BibbyPet bottle top dispenser provides a low cost option for liquid dispensing without compromising on accuracy or reliability. The Bibbypet's innovative design utilises materials suitable for routinely dispensing many common laboratory liquids; dilute acids, buffers, biological solutions, polar solvents, even strong alkaline solutions.


\* The values of accuracy and coefficient of variation are final test values referring to the delivered nominal volume, instrument and distilled water at equilibrium with ambient temperature (20°C) and with smooth operation.

| Model | Volume, ml | Subdivision, ml | Accuracy, % | Coefficient of variation, % |
|-------|------------|-----------------|-------------|-----------------------------|
| BD10  | 1 to 10    | 0.2             | 1.2         | 0.2                         |
| BD25  | 2.5 to 25  | 0.5             | 1.2         | 0.2                         |

### Accessories


#### Polypropylene adapters

| Model  | Adapter thread, mm | For bottle thread, mm | Pack Qty |
|--------|--------------------|-----------------------|----------|
| BD3222 | GL32               | GL22                  | 1        |
| BD3228 | GL32               | GL28                  | 1        |
| BD3230 | GL32               | GL30                  | 1        |
| BD4532 | GL45               | GL32                  | 1        |
| BD4538 | GL45               | GL38                  | 1        |
| BD3219 | GL32               | 19/26                 | 1        |
| BD3224 | GL32               | 24/29                 | 1        |
| BD3229 | GL32               | 29/32                 | 1        |



# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)



Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)

## Dispensing cartridge

| Model  | Capacity | Materials              | Pack Qty |
|--------|----------|------------------------|----------|
| BD1010 | 10ml     | Piston PE, cylinder PP | 3        |
| BD1012 | 25ml     | Piston PE, cylinder PP | 3        |

## Replacement parts

| Model  | Description         | Pack Qty |
|--------|---------------------|----------|
| BD1014 | Discharge tube      | 1        |
| BD1016 | Filling tube, 250mm | 2        |
| BD1018 | Filling tube, 500mm | 2        |
| BD1026 | Seals               | 6        |



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



see the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

