

SA8



Mixer, vortex, variable speed,

Mixing speed can be selected from 200rpm for very gentle mixing to 2500rpm for vigorous agitation. Choose between 'touch' mode or continuous operation simply by pressing a button on the fascia. Robust die-cast body avoids unnecessary movement during use. The ergonomic low profile design makes everyday vortexing comfortable for the user.

Integral retort rod fixing allows vessels to be secured above the vortex action for long-term mixing. The SA8 can mix a variety of other vessels when used in conjunction with the SA8/1 accessory pack detailed below.

Technical Specification

| | |
|----------------------------|---------------------------|
| Speed | 200 to 2500rpm |
| Orbit diameter, mm | 4.2 |
| Touch mode | Yes |
| Continuous mode | Yes |
| Dimensions, mm (w x d x h) | 135 x 215 x 78 |
| Net weight, kg | 3.2 |
| Electrical supply | 90 – 240V, 50 / 60Hz, 20W |
| IP Rating | 31 |

Key Features

- Variable speed control between 200 and 2500rpm
- Intermittent or continuous mode
- Stable low profile body
- Robust die-cast construction



SA8

Ordering Information

| Model | Description |
|-------|------------------------------|
| SA8 | Vortex mixer, variable speed |

Accessory pack

For use with SA8 above, pack contains:

- 1 x Plastic cradle for standard microtitre plate
- 1 x Circular foam insert for beakers and flasks (up to 500ml capacity)
- 1 x Rectangular foam insert (holds 8 x 0.2ml, 8 x 0.5ml and 16 x 1.5ml microcentrifuge tubes)
- 1 x Solid rectangular foam insert (blank for custom holes)

Ordering Information

| Model | Description |
|-------|--|
| SA8/1 | Accessory pack (for use with SA8 only) |



SA8 with cradle

with SA8/1 Microtitre plate cradle (plate not included)

Key Features

- Fixed speed 2500rpm
- Robust die-cast construction
- Automatic press start



SA7



Mixer, vortex, fixed speed,

For rapid mixing of samples contained in test-tubes, small flasks and bottles. Starts automatically when the rubber cup is depressed and stops once vessel is removed.

Heavy die-cast body for stability during use.

Technical Specification

| | |
|----------------------------|---------------------------|
| Speed | Fixed 2500rpm |
| Orbit diameter, mm | 4.2 |
| Dimensions, mm (w x d x h) | 135 x 215 x 78 |
| Net weight, kg | 3.2 |
| Electrical supply | 90 – 240V, 50 / 60Hz, 20W |
| IP Rating | 31 |

Ordering Information

| Model | Description |
|-------|---------------------------|
| SA7 | Vortex mixer, fixed speed |



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax:01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

