



## Calibrate with confidence using the Tecal “Field” Range of Dri-Block Calibrators

A new range of lightweight dry-block thermal calibrators that operate from 45°C below ambient up to 650°C has been developed by Techne for convenient use in the field.

The three models in the new Tecal F range incorporate the precision control system used on Techne’s Tecal H and S models to provide excellent reliability and accuracy, with a number of improvements such as an enhanced display resolution of 0.01°C and faster heating and cooling rates.

The Tecal 140F operates between 45°C below ambient (typically -20°C) and 140°C with a temperature accuracy of  $\pm 0.3^\circ\text{C}$  and a temperature uniformity of  $\pm 0.2^\circ\text{C}$ . Heating time to 100°C is five minutes.

All the models are offered with or without RS-232 communications and can be used with a wide range of inserts to suit different probe sizes. Automatic fan cooling is standard, with accelerated cooling probes also available.

## Tecal 140F

### Technical Specification

Low temperature Portable Field Temperature Calibrator	
Minimum temperature	45°C below ambient (typically -20°C in ambient of 25°C)
Maximum temperature	140°C
Temperature accuracy	$\pm 0.3^\circ\text{C}$
Temperature uniformity	$\pm 0.2^\circ\text{C}$
Temperature stability	$\pm 0.05^\circ\text{C}$ (after 10 mins)
Display resolution	0.01°C or 0.1°F
Set point resolution	0.1°C or 0.1°F
Heating rate, 20°C to 100°C	5 minutes
Cooling rate, 100°C to 0°C	9 minutes
Large well	Ø38mm x 114 mm inserts
Fan cooling	Automatic
Dimensions HxWxDmm	273 x 207 x 289
Weight	11kg

 See page 9 for full range of inserts available

### Options

- Inserts from standard list type FINSAL-
- RS-232 interface plus software and cable
- Carry case soft
- Carry case hard
- UKAS calibration certificate

### Ordering Information

Product Code	Model	Voltage	Hz	Watts	Net Weight (kg)
FDB140FD	Tecal 140F	230v	50/60	400	11
FDB140FP	Tecal 140F	120v	50/60	400	11
FDB140FY	Tecal 140F	100v	50/60	400	11
FDB140FR	Tecal 140F with RS-232	230v	50/60	400	11
FDB140FS	Tecal 140F with RS-232	120v	50/60	400	11
FDB140FT	Tecal 140F with RS-232	100v	50/60	400	11
<b>Accessories</b>					
FCALFCSL	Hard Carry case for 140F				8.2
FSC140F	Soft Carry case for 140 + 650F				2.3
6104266	RS-232 cable for use relevant units				0.2



140F



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

