

# WOLF LABORATORIES

Tel: 01759 301142 Fax: 01759 301143 Email: [sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)  
Website: [www.wolflabs.co.uk](http://www.wolflabs.co.uk)

## Cryo Boxes



Cardboard boxes WxD 133x133mm

All boxes listed below have an external measure of 133x133 mm. Our boxes are offered in 3 different heights, 50, 75 and 100 mm.

The boxes are made of fibreboard/cardboard and are Water-repellent for longer durability.

The table below shows the article number by height.

Height	Article no. Box 50 mm	Article no. Box 75 mm	Article no. Box 100 mm
White	TE22005	TE22050	TE22100



Grid dividers

For the boxes above we offer 13 different grid dividers. Choose the number of cells and for which box height the divider should be used.

The 30 mm is for 50 mm boxes  
The 40 mm is for 75 mm boxes  
The 65 mm is for 100 mm boxes

Grid dividers, height	Article no. 30 mm for 50 mm box	Article no. 40 mm for 75 mm box	Article no. 65 mm for 100 mm box
5x5 cells for vials up to 25 mm O.D.	TE22111	-	-
7x7 cells for vials up to 16 mm O.D.	TE22106	TE22130	TE22140
8x8 cells for vials up to 14 mm O.D.	TE22108	TE22132	TE22142
9x9 cells for vials up to 13 mm O.D.	TE22110	ZXZTE22134	TE22144
10x10 cells for vials up to 12 mm O.D.	TE22120	TE22136	TE22146



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.