

WOLF LABORATORIES

Tel: 01759 301142 Fax: 01759 301143 Email: sales@wolflabs.co.uk
Website: www.wolflabs.co.uk

Cryo Boxes



Cardboard boxes WxD 133x133mm

All boxes listed below have an external measure of 133x133 mm. Our boxes are offered in 3 different heights, 50, 75 and 100 mm.

The boxes are made of fibreboard/cardboard and are Water-repellent for longer durability.

The table below shows the article number by height.

| Height | Article no. Box 50 mm | Article no. Box 75 mm | Article no. Box 100 mm |
|--------|--------------------------|--------------------------|---------------------------|
| White | TE22005 | TE22050 | TE22100 |



Grid dividers

For the boxes above we offer 13 different grid dividers. Choose the number of cells and for which box height the divider should be used.

The 30 mm is for 50 mm boxes
The 40 mm is for 75 mm boxes
The 65 mm is for 100 mm boxes

| Grid dividers, height | Article no. 30 mm for 50 mm box | Article no. 40 mm for 75 mm box | Article no. 65 mm for 100 mm box |
|--|------------------------------------|------------------------------------|-------------------------------------|
| 5x5 cells for vials up to 25 mm O.D. | TE22111 | - | - |
| 7x7 cells for vials up to 16 mm O.D. | TE22106 | TE22130 | TE22140 |
| 8x8 cells for vials up to 14 mm O.D. | TE22108 | TE22132 | TE22142 |
| 9x9 cells for vials up to 13 mm O.D. | TE22110 | ZXZTE22134 | TE22144 |
| 10x10 cells for vials up to 12 mm O.D. | TE22120 | TE22136 | TE22146 |



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

