

Explosion-Proof Lab Refrigerators, Freezers and Combination Units

Designed to be used in areas where ignitable materials are normally present.

Thermo Scientific explosion-proof laboratory refrigerators, freezers and combinations units are particularly suited for safe storage in hazardous environments.

Not available in North America, consult your sales representative for availability.

- Sparkfree interior and exteriors reduce the risk of explosion
- All models are manual defrost
- Tough white exterior
- CFC- and HCFC-free refrigeration system and insulation
- Designed for use in CLASS I, GROUPS C & D & CLASS I, ZONE1, GROUP IIB
- Plug-free cords (units must be hardwired to wall)

Warranty: 13 months parts and labor



Model No: 20EFEETSV
591.82 liter explosion-proof freezer



Model No: 10ECEETSV
286 liter explosion-proof
combination refrigerator/freezer

Thermo Scientific Explosion-Proof Refrigerators (Temperature Range 1° to 12°C)

| Model No. | Capacity | Interior Storage | Voltage | Exterior Dimensions: (H x W x D) | Interior Dimensions: (H x W x D) | Shipping Weight |
|------------|---------------------------------|------------------------------------|-----------|--|--|-----------------------|
| 20ERCETSV* | 20.9 cu. ft. (591.82 liters) | 4 fixed shelves | 230V/50hz | 77.5" x 34" x 34.96" (196.9 x 86.3 x 88.8 cm) | 57.8" x 27.8" x 23.8" (146.7 x 70.49 x 60.3 cm) | 198 lbs. (89.8 kg) |
| 20EREETSV | 20.9 cu. ft. (591.82 liters) | 4 fixed shelves, 6 door shelves | 230V/50hz | 77.5" x 34" x 34.96" (196.9 x 86.3 x 88.8 cm) | 57.8" x 27.8" x 23.8" (146.7 x 70.49 x 60.3 cm) | 198 lbs. (89.8 kg) |

* Corrosion-resistant model, features stainless steel interior door and stainless steel shelves

Thermo Scientific Explosion-Proof Freezers (Temperature Range -12° to -20°C)

| Model No. | Capacity | Interior Storage | Voltage | Exterior Dimensions: (H x W x D) | Interior Dimensions: (H x W x D) | Shipping Weight |
|-----------|---------------------------------|------------------------------------|-----------|--|--|-----------------------|
| 20EFEETSV | 20.9 cu. ft. (591.82 liters) | 4 fixed shelves, 6 door shelves | 230V/50hz | 77.5" x 34" x 34.96" (196.9 x 86.3 x 88.8 cm) | 57.8" x 27.8" x 23.8" (146.7 x 70.49 x 60.3 cm) | 198 lbs. (89.8 kg) |

Thermo Scientific Explosion-Proof Combination Refrigerator/Freezers (Temperature Range: 1° to 12°C [ref], -12° to -20°C [frz])

| Model No. | Capacity | Interior Storage | Voltage | Exterior Dimensions: (H x W x D) | Interior Dimensions: (H x W x D) | Shipping Weight |
|-----------|--------------------------------|---|-----------|--|--|---------------------|
| 10ECEETSV | 10.1 cu. ft. (286.0 liters) | 3 adjustable interior, 1 fixed interior, 1 drawer, 4 door shelves | 230V/50hz | 59.75" x 23.63" x 28.6" (151.77 x 60 x 76.2 cm) | Ref: 38.25" x 19.75" x 18.25" (97.2 x 50.2 x 46.4 cm); Frz: 11.25" x 18" x 18.25" (28.6 x 45.7 x 46.4 cm) | 150 lbs. (68 kg) |

ThermoFisher
SCIENTIFIC



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.