

U100H

1.5 Litre Ultrasonic Bath (Heated)

Ultrawave's **U-Series** is compact and user friendly, offering effective cleaning and processing of small items in both professional and domestic applications.

These include healthcare, podiatry, tattooists, jewellers and many more.

Available in heated and unheated versions, the **U-Series** allows simple control of time (0 to 15 minutes) and temperature (ambient to 70°C) where heated options are selected.



The **U100H** is a 1.5 Litre ultrasonic bath with heating function to reduce processing time.

Each **U-Series** is supplied complete with a stainless steel basket and ABS plastic lid, useable as a drip tray to collect excess water when the basket is removed from the tank.

The **U100H** provides the cost benefits of excellent entry level ultrasonics with Ultrawave's long renowned quality and reliability.



U100H

1.5Litre Ultrasonic Bath (Heated)

Construction Case	304 Stainless Steel
Construction Tank	304 Stainless Steel
Control	Analogue (Time)
Overall Dimensions	180 x 160 x 190mm (WxDxH)
Tank Dimensions	150 x 140 x 100mm (WxDxH)
Basket Dimensions	120 x 95 x 55mm (WxDxH) (75mm height with handles)
Maximum Capacity	1.75Litres
Working Capacity	1.5Litres
Ultrasonic Power	35W
Ultrasonic Power (per litre)	23.3W/litre
Operating Frequency	44KHz
Maximum Heating Capacity	Ambient - 70°C *
Heater Power	80W
Timer	0 - 15 minutes
Mains Connection	220 – 240V (50 - 60Hz)
Drain Outlet	N/A
Lid	ABS Plastic
Basket	Stainless Steel
Gross Weight (inc. lid, basket & packaging)	4.5KGs
Net Weight	3.7KGs
Warranty	3 Years

*Temperature is non-adjustable



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

