

U50

0.5Litre Ultrasonic Bath

Ultrawave's **U-Series** is compact and user friendly, offering effective cleaning and processing of small items in both professional and domestic applications.

These include healthcare, podiatry, tattooists, jewellers and many more.

Available in heated and unheated versions, the **U-Series** allows simple control of time (0 to 15 minutes) and temperature (ambient to 70°C) where heated options are selected.



The **U50** is a 0.5Litre ultrasonic bath with a timer control of up to 15 minutes.

Each **U-Series** is supplied complete with a stainless steel basket and ABS plastic lid, useable as a drip tray to collect excess water when the basket is removed from the tank.

The **U50** provides the cost benefits of excellent entry level ultrasonics with Ultrawave's long renowned quality and reliability.



U50

0.5Litre Ultrasonic Bath

| | |
|--|---|
| Construction Case | 304 Stainless Steel |
| Construction Tank | 304 Stainless Steel |
| Control | Analogue (Time) |
| Overall Dimensions | 160 x 95 x 165mm (WxDxH) |
| Tank Dimensions | 150 x 85 x 65mm (WxDxH) |
| Basket Dimensions | 120 x 60 x 25mm (WxDxH) (50mm height with handles) |
| Maximum Capacity | 0.6Litres |
| Working Capacity | 0.5Litres |
| Ultrasonic Power | 35W |
| Ultrasonic Power (per litre) | 70W/litre |
| Operating Frequency | 44KHz |
| Maximum Heating Capacity | N/A |
| Heater Power | N/A |
| Timer | 0 - 15 minutes |
| Mains Connection | 220 – 240V (50 - 60Hz) |
| Drain Outlet | N/A |
| Lid | ABS Plastic |
| Basket | Stainless Steel |
| Gross Weight (inc. lid, basket & packaging) | 1.8KGs |
| Net Weight | 1.4KGs |
| Warranty | 3 Years |



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

