



# DOC-PRINT EASY & FUN

*Stand-Alone Gel Documentation Imaging*

## Camera & Optics

CX3 camera:

- Scientific grade camera - Made in Germany
- 3 megapixels resolution extendable to 20 megapixels
- Passive cooling
- 16-bit – 65 536 grey levels
- USB-3 connection
- Field of view: 24x18cm

## Ease Of Use

- One-Click-to-the-Image™
- Auto-exposure
- Self-understandable user interface
- Print or save on an external drive via the USB port
- Connect to a thermal printer to get vibrant glossy printouts

## Hardware

ABS hood  
F590 emission filter for standard gel applications  
Transilluminator not included.  
Complete network connectivity

- Very large 12.6" touch screen interface display with Magster™ technology
- Stand-alone system - No computer require
- Small footprint

## Software

Free VisionCapt software with full GLP compliance.  
Molecular weight calculation, band quantification, colony counting, distance calculation, text annotation and image enhancement included.

## Applications

DNA and RNA gels and fluorescence stain imaging:

- Ethidium bromide, Sybr-Safe, Sybr-Green, Gel-Red, Gel-Green, Sybr-Gold, GFP, Pro-Q Emerald, Sypro ruby, FITC, DAPI

Colorimetric stained protein gels, X-Ray film, autorads, SSCP gels, colony dish and flask imaging and other EPI white light applications:

- Coomassie blue, Silver stain, Ponceau S Red, Copper stain...

Sky- Pad or blue light conversion screen for DNA/RNA detection (avoid « nicking » DNA):

- Sybr Safe, Sybr Green, eGFP

## Technologies & Innovation

3D Dynamics Scan  
Magster™ display technology

## Options

Several different Application-Pad can be easily used with the Doc-Print:

- UV-Pad – 312nm or dual wavelength 312 and 365nm
- Super-Bright-Pad – 312nm or dual wavelength 312 and 365nm
- Sky-Pad: LED blue light transilluminator, 470nm
- White-Light- Pad: white light LED transilluminator

White light or blue light conversion screen  
Advanced Bio-1D quantification software



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.