

SIDE ACCESS RACKS

Side Opening Horizontal units also known as side access racks, fit onto shelves of upright freezers to hold cryoboxes or microtiter plates. Side access racks feature a handle on one end to remove it from the freezer, these racks are best suited for longer term storage. Side access racks are manufactured from aluminium or stainless steel. Side access racks can be any size to meet the requirements of the user. They can also feature a locking rod to hold items in place, and have a handle fitted to the top to make the unit easier to transport.

Dimensions (mm)	Box Size / Boxes Per Rack	Part Number
295H x 145W x 290D	4" (100mm) / 4	SOR455
295H x 145W x 290D	3" (75mm) / 6	SOR453
295H x 145W x 290D	2" (50mm) / 10	SOR110
420H x 140W x 560D	4" (100mm) / 12	SOR304
420H x 140W x 560D	3" (75mm) / 20	SOR368
420H x 140W x 560D	2" (50mm) / 28	SOR090
235H x 141W x 560D	2" (50mm) / 16	SOR196
235H x 141W x 560D	3" (75mm) / 12	SOR197
290H x 140W x 560D	2" (50mm) / 20	SOR032
340H x 140W x 560D	2" (50mm) / 24	SOR033
230H x 140W x 560D	2" (50mm) / 16	SOR196

Also available in stainless steel.


Wesbart also offers a very popular adjustable side access rack which offers the user flexibility to store different size boxes. Using the unique design of this rack the user can set the removable shelves to store either 16x 2" boxes, 12x 3" boxes or 8x 4" boxes. *Part Number SOR002*

BESPOKE RACKS


Should you require a bespoke front or side access rack please visit our website (www.wesbart.co.uk) to fill in a quotation request. The information that will be required is the freezer make and model, maximum internal size of the shelf and the size of the item you require racking.



Side access rack to store 2" box **SOR196**



Side access rack to store microtiter plates **SOR004**



Side access rack to store various height boxes **SOR002B**

EPPENDORF F740 (3 SHELF)

Wesbart offers a large range of racking to make the most of the space available with the Eppendorf F740. This new freezer will hold up to 576 standard 2inch cryoboxes. That's 57,600 2ml samples 7,200 more than the U725.

Below are details of the racking to store 2, 3 or 4 inch cryoboxes along with details for horizontal drawer systems and storage bins, ideal for organising loose items. Wesbart racks are available in aluminium or stainless steel.

We pride ourselves on the quality of our workmanship and products. All our racks have no sharp edges or corners and are designed to last the lifetime of your freezer. We are so confident with the quality and durability of our products we are the only racking manufacturer to offer a 3-Year warranty.

RACKS FOR CRYOBXES

Box Size	Boxes Per Rack	Racks Per Shelf	Racks Per Freezer	Part Number (Front Access Rack)	Part Number (Side Access Rack)
2" (50mm)	32	6	18	BH248/8T	SOR500
3" (75mm)	20	6	18	AH179/5UT	SOR368
4" (100mm)	16	6	18	BH248/4UT	SOR501

STORAGE BINS

Bins Per Shelf	Bins Per Freezer	Dimensions (mm)	Part Number
4	12	410H x 205W x 560D	BIN 595

HORIZONTAL DRAWER SYSTEMS

Units Per Shelf	Drawer Per Rack	Dimensions (mm)	Part Number
6	2	420H x 140W x 560D	AH179/2DR
6	3	420H x 140W x 560D	AH179/3DR
4	2	420H x 205W x 560D	AH178/2DR
4	3	420H x 205W x 560D	AH178/3DR

All racking is suitable for cryoboxes with an external footprint of 133 x 133mm. All front access racks include trays with handgrips.

If you require bespoke racking for the storage of centrifuge tubes, bottle and microtiter plates please contact us for a quotation. Whatever your racking requirements we will have a solution for you to maximise the use of your freezer.



AH179/5UT




BH248/8T




AH178/2DR



BH249/6T



SOR196




AH506/4T

EPPENDORF F740 (5 SHELF)

Wesbart offers a large range of racking to make the most of the space available with the Eppendorf F740. With the 5-shelf option, you are still able to store 528 standard 2 inch cryoboxes, 24 more than the U725. With the shelf arrangement of this freezer it is best suited to store a range of different height cryoboxes. With the bottom shelf (shelf 5) holding a different size rack.

Below are details on racking to store 2, 3 or 4 inch cryoboxes along with details for horizontal drawer systems and storage bins, ideal for organising loose items. Wesbart racks are available in aluminium or stainless steel.

We pride ourselves on the quality of our workmanship and products. All our racks have no sharp edges or corners and are designed to last the lifetime of your freezer. We are so confident with the quality and durability of our products we are the only racking manufacturer to offer a 3-Year warranty.



RACKS FOR CRYOBXES

SHELF 1, 2, 3, 4

Box Size	Boxes Per Rack	Racks Per Shelf	Racks Per Freezer	Part Number (Front Access Rack)	Part Number (Side Access Rack)
2" (50mm)	16	6	24	AH506/4T	SOR196
3" (75mm)	12	6	24	AH506/3T	SOR197
4" (100mm)	8	6	24	AH506/2T	SOR198

SHELF 5

Box Size	Boxes Per Rack	Racks Per Shelf	Racks Per Freezer	Part Number (Front Access Rack)	Part Number (Side Access Rack)
2" (50mm)	24	6	6	BH249/6T	SOR502
3" (75mm)	16	6	6	BH249/4T	SOR503
4" (100mm)	12	6	6	BH249/3T	SOR504

STORAGE BINS

Bins Per Shelf	Bins Per Freezer	Dimensions (mm)	Part Number
6	30	235H x 140W x 560D	BIN853
6	30	235H x 140W x 560D (with lid)	BIN853L

HORIZONTAL DRAWER SYSTEMS

Units Per Shelf	Drawer Per Rack	Dimensions (mm)	Part Number
3	2	235H x 280W x 560D	GH90/2DR

All racking is suitable for cryoboxes with an external footprint of 133 x 133mm. All front access racks include trays with handgrips.

If you require bespoke racking for the storage of centrifuge tubes, bottle and microtiter plates please contact us for a quotation. Whatever your racking requirements we will have a solution for you to maximise the use of your freezer.



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

