

Whatman Membrane-Butler



The Membrane-Butler offers optimal handling of Whatman ME and MicroPlus CN STL membranes. With each turn a membrane filter is ejected from its sterile packaging and can be removed easily with tweezers.

- Durable design for reliable performance
- Easy grasp handle to dispense filter
- Minimised cross-contamination by eliminating handling of individually sealed filters
- Suitable for use on sterile benches
- Compact size for portability and preserve area of workspace



Code	Alt Ref	Description	Price
FIL5475	10477100	Manual Membrane-Butler	£1116.00

Corning® PES Vacuum Filtration Systems



- Corning's vacuum filter products use a square filter membrane support system that increases filtration rates
- Improved, higher capacity polyethersulfone (PES) membranes maximise filter throughput

Vacuum Filter Systems

- Low liquid retention funnel design
- Filters feature printing on the funnel for easy product identification
- Angled hose connector simplifies vacuum line attachment
- Receiver bottles feature easy grip sides for improved handling
- Extra plastic storage bottles are available
- Prefilters must be ordered separately
- Sterile and nonpyrogenic
- PES Membrane
- Caps for receiver bottles are sterile and individually packaged



Supplied in packs of 12

Code	Capacity, mL	Pore Size, µm	Membrane, cm ²	Price
431153	150	0.22	13.6	£98.00
431096	250	0.22	33.2	£101.00
431097	500	0.22	33.2	£166.00
431098	1000	0.22	54.5	£241.00

Bottle Top Vacuum Filters

- Low liquid retention funnel design
- Available in 33mm and 45mm neck sizes to fit most glass and plastic media storage bottles
- 45mm neck sizes fit on Corning® glass 1395 bottles, to be used on 2L and smaller sizes only
- Prefilters must be ordered separately
- Sterile and nonpyrogenic
- PES membrane



Code	Capacity, mL	Fits Bottle Neck Size, mm	Pore Size, µm	Membrane, cm ²	Pack	Price
431160	150	33	0.22	13.6	48	£284.00
431161	150	45	0.22	13.6	48	£284.00
431117	500	33	0.22	33.2	12	£106.00
431118	500	45	0.22	33.2	12	£106.00
431174	1000	45	0.22	54.5	12	£178.00



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

