

**REFERENCE :**

**062646CL**

**Rotating mixer ClearLine®**

**Features :**

- Rotating shaker for tubes from 1.5 to 50 ml
- Adjustable stirring speed and angle
- Reversible mode: 360° rotation then change of direction
- Comes with 3 carousels :
  - for 48 micro tubes of 1.5 or 2.0 ml
  - for 20 tubes of 15 ml
  - for 12 tubes of 50 ml



Speed (rpm)	1 to 40
Angle	0 to 90°
Timer	1sec to 9h
Capacity	48 microtubes de 1,5ml 20 tubes de 15 ml 12 tubes de 50 ml
Dimensions L x P x H (mm)	210 x 270 x 320
Temperature	4 to 60°C
Weight (kg)	3,2
Power supply	230V, 50-60Hz

**Unité de vente :**

1 agitateur



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.