

## Tilt Shaker WS 10



### Description

The Tilt Shaker WS 10 works with slow, linear tilting motion and a fix angle of tilt of 5° with a shaking speed of 1 - 23 rpm. The shaking platform of 340 x 340 mm can be loaded with a weight of 10 kg.

### Basic Equipment

Basic device incl. basic tray and rubber mat

### Details

**Motion:** linear tilting motion

**Shaking platform:** 340 x 340 mm

**Max. load:** 10 kg

**Shaking speed:** 1 - 23 rpm, steplessly adjustable  
**Angle of tilt:** 5° fixed  
**Timer:** 0 - 120 min / continuous  
**Order No.:** 6114 000  
**Accessory:** Multi-storey rack system WS / TL

## Technical specifications

**Electrical supply:** 24 V DC, 50 / 60 Hz 115 V version on request  
**Enclosure protection:** IP 21  
**Heat emission:** approx. 8 W  
**Ambient temperature:** 5°C to 50°C  
**Relative humidity:** ~ 85 %  
**Dimensions (w x d x h):** 355 x 450 x 205 mm  
**Weight:** 18,5 kg



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.