

# Compact Shaker KS 15 A control



## Description

The programmable Orbital Shaker KS 15 A control with digital display for reproducible shaking tasks. The loading capacity is 15 kg on a shaking platform of 400 x 300 mm.

### Basic Equipment

Basic device without shaking plate, **without rack system**

With digital display and programmable control.

### Loading Capacity

Erlenmeyer flasks 100 ml	20 pieces
Erlenmeyer flasks 250 ml	12 pieces
Erlenmeyer flasks 500 ml	6 pieces
Erlenmeyer flasks 1000 ml	4 pieces
Separating funnels 250 ml	1 pieces

You will find further loading possibilities [here](#)

**Option:**[Control interface](#)[USB interface](#)

## Details

<b>Motion:</b>	orbital
<b>Shaking platform:</b>	400 x 300 mm
<b>Max. load:</b>	15 kg
<b>Shaking speed:</b>	30 - 420 rpm, in steps of 5
<b>Stroke:</b>	17 mm
<b>Runtime:</b>	programmable / continuous
<b>Order No.:</b>	6170 000
<b>Accessory:</b>	Rack System Combifix KS Universal tray KS Incubator Hood TH 15

## Technical specifications

<b>Electrical supply:</b>	230 V or 115 V, 50 / 60 H please indicate in case of order
<b>Enclosure protection:</b>	IP 21
<b>Heat emission:</b>	approx. 10 - 15 W
<b>Ambient temperature:</b>	5°C to 50°C
<b>Relative humidity:</b>	~ 85 %
<b>Dimensions (w x d x h):</b>	510 x 490 x 150 mm
<b>Weight:</b>	19 kg

---



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.