



LABCOLD™

medical & scientific refrigeration



RLVF1517-ATEX Free Standing Sparkfree ATEX Freezer

Part Number	RLVF1517- ATEX
Capacity (Litres)	406
Style	Free standing
Exterior Dimensions (HxWxD)	1875 x 600 x 700mm
Interior Dimensions (HxWxD)	1690 x 440 x 500mm
Net Weight	97kg
Energy Consumption	3.16kWh/24hr
Operating Temperature	-18°C to -25°C
Factory Set Point	-20°C
Interior Construction	Painted Aluminium
Exterior Construction	Painted Steel
Door Style	Solid
Door Hinging	Right
Door Earth Connection	✓
Door Lock	✓
Door Open Alarm	✓
Adjustable Shelves	6
Shelf dimensions (WxD)	413mm x 438mm
Fan Circulated (ATEX fan)	✓
Rollers	4
Temperature Display	✓
High Temperature Alarm	✓
Low Temperature Alarm	✓
Automatic Defrost	✗
2 Year Warranty (UK only)	✓

The Labcold ATEX freezer range has been designed with an interior which is licensed as a Zone hazardous area under 2014/34/EU (ATEX Directive) in order to store products which are classed in explosion group IIB+H2 and temperature class T6 under 2014/34/EU (ATEX Directive).

Warning: Do not store open containers of volatile substances in this freezer.

CLAIM
YOUR
FREE
YEAR
WARRANTY



All Labcold sparkfree ATEX freezers are fitted with a cylinder lock with two keys and a controller with a digital temperature display. The controller also features a high/low audio and visual temperature alarm in addition to a door open alarm.





WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.