

## Heating Stirrer Mantle (20,000ml)

**Product Code:** ES-309N

**Brand:** Medline Scientific

This mantle features an integral magnetic stirrer, which is independently controlled. The stirrer speed may be varied from 100 to 1000rpm and the temperature is regulated by an infinity bimetallic thermostat. Retort rod holder integrated into base of unit to accommodate our range of accessories.

### Specifications:

Volume (ml)	20,000
Case Construction	Steel, with chemical resistant coating
Thermal Insulation	Moulded silica powder
Heating Element	Nichrome wire
Mantle	Woven glass fibre mat
Element Temperature	Up to 450 °C
Speed Range	10 - 1000rpm
Power Supply	220V AC
Watts	1800





# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.