

## Medline Analogue Hotplate (MS-100)

**Product Code:** MS-100

**Brand:** Medline Scientific

Medline's Analogue Hotplate (MS-100) offers a large heating area of 180 x 180mm with safe controlled heating up to 400°C. Constructed from chemically resistant ceramic coated steel, this lab hotplate delivers excellent reliability with reproducible results at a sensible price. A standard 12.5mm screwed thread is included for the addition of a retort rod (not included).



## Specifications

|                        |                   |
|------------------------|-------------------|
| Description            | Hotplate          |
| Plate Area (mm)        | 180 x 180         |
| Temperature Range (°C) | Max. 400          |
| Dimensions (W x D x H) | 410 x 210 x 150mm |
| Weight (kg)            | 4                 |



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.