


CAT. NO.....	295-3345		295-4545	
PORE SIZE (µm).....	0.45		0.45	
UPPER CAP. (mL).....	500		500	
BOTTLE NECK SIZE....	33mm		45mm	
MEMBRANE COLOR.....	WHITE		WHITE	
MEMBRANE DIA (mm)...	75		75	
	INCHES	MM	INCHES	MM
A (UPPER I.D.).....	3.8	96	3.8	96
B HEIGHT.....	5.0	127	5.1	130
C DIA (BASE).....	1.4	36	1.9	48
D DIM.....	2.5	64	3.0	76
NOM. WT. (gms).....	70		72	

MATERIALS	
UPPER CONTAINER	PS
COVER	PS
COLLAR	BLUE PS
ADAPTOR	WHITE LOPE
GASKET	FOAMED LOPE
MEMBRANE	PES
VENT PLUG	CA



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
**SCIENTIFIC**

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	MC <sup>2</sup>	APP'D	PVC	PSS 00-0295-XX45	REV. 1205
CK'D	SH		12/22/05		
DESCRIPTION:					
BOTTLE TOP FILTERS					



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.