

# Chest freezer XLT C150



- Low cost of ownership
- Limited maintenance and on-site service
- Easily installable & movable
- Efficient cyclopentane insulation
- Eco-friendly, HCFC-free Nature R refrigerant
- ISO 9001, 14001 & 18001 OHSAS Certified

## Specifications

### Technical

|                          |                  |
|--------------------------|------------------|
| Gross volume             | 140 litres       |
| Net Volume               | 133 litres       |
| Temp. range              | -45°C to -65°C   |
| Temp. probes             | 1/PT 1000        |
| Compressor               | Embraco          |
| Baskets (incl./optional) | 1/2              |
| Refrigerant              | Nature R         |
| Noise level              | 55 dB            |
| Power consumption 24/h   | 3,7 kWh          |
| Power consumption        | 1.350 kWh / year |
| Voltage                  | 230V             |
| Frequency                | 50 Hz            |

### Features

|  |        |
|--|--------|
| Audible and visual alarm                 | YES    |
| Castors                                  | 4      |
| Integrated lock                          | YES    |
| Inner lids                               | NO     |
| Defrost                                  | Manual |
| Battery Backup                           | NO     |
| USB                                      | NO     |
| Data connection                          | NO     |
| Digital controller with display          | XR30CX |
| Data logging                             | NO     |
| High/Low temperature alarm               | YES    |
| Porthole for external probe & CO2 backup | Ø12 mm |
| Prepared for external logging equipment  | YES    |

### Material

|       |                             |
|-------|-----------------------------|
| Inner | Precoated white steel       |
| Outer | Galvanized pre-coated steel |

### Outer dimensions

|                                  |        |
|----------------------------------|--------|
| Height (with castors)            | 890 mm |
| Width                            | 720 mm |
| Depth (incl. condenser)          | 630 mm |
| Depth (incl. condenser & handle) | 680 mm |
| Net weight                       | 56 kg  |

### Inner dimensions

|                   |            |
|-------------------|------------|
| Height            | 635/385 mm |
| Width             | 350/210 mm |
| Depth             | 440 mm     |
| Insulation        | 80 mm      |
| Insulation type   | PUR        |
| Internal humidity | <5 %       |

## XLT C150



Product no.

N112002

## Storage capacity

| Box height | No. of boxes | No. of racks | Total no. of boxes | Product no. |
|------------|--------------|--------------|--------------------|-------------|
| 50 mm      | 6            | 3            | 18                 | N311200     |
|            | 11           | 6            | 66                 | N311220     |
|            |              |              | <b>= 84</b>        |             |
| 75 mm      | 4            | 3            | 12                 | N311300     |
|            | 8            | 6            | 48                 | N311320     |
|            |              |              | <b>= 60</b>        |             |
| 100 mm     | 3            | 3            | 9                  | N311400     |
|            | 5            | 6            | 30                 | N311420     |
|            |              |              | <b>= 39</b>        |             |



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.