



WolfLabs

PURSUIT

2026 Issue **02**

July - December



Haier
Class II
Cabinets

38%
Off

Grant ORCA
Recirculating chiller



25%
Off

Hettich
Centrifuges



45%
Off

Haier
UltraECO ULT Freezers



46%
Off

Jump to

Index

Aspirators	44	Microbalances	19
Autoclaves	9-10	Microscopes	39-41
Balances - Analytical	17-19	Ovens - Drying/Warming	27
Balances - Portable	21	Ovens - Vacuum	28
Balances - Top Pan	20	pH and Conductivity Meters	37
Block Heaters	45	Pipettes - Manual	35-36
Bunsen Burners	38	Pipetting Controllers	36
Centrifuges and Microfuges	5-6	Pumps - Syringe	29
Chairs - Laboratory	34	Pumps - Vacuum - General Purpose	28
Chillers - Recirculating and Water Baths - Cooled	25-26	Rockers	43
Class II Cabinets	3	Rollers	43
CO2 Incubators	7	Rotators	43
Cooled Incubators	8	Servicing	16
Cryostats	48	Shakers - Microplate	42
Cryostorage Equipment - Liquid Nitrogen Sample Storage	4	Shaking Incubators	9
Density Meters	38	Sonicators - Probe Type	34
Freezers - Ultra-Low Temperature -80°C	2-3	Spectrophotometers	29
Freezers -20°C	11	Stirrers - Hotplate	32
Fridge Freezers	13	Stirrers - Magnetic	31-32
Fridges	12	Stirrers - Overhead	30
Fridges - Pharmacy	14	Storage Cabinets	42
Fume Cupboards - Recirculating	23	Thermal Cyclers	46
Gel Documentation Systems	38	Thermoshakers	47
Glass Washers - Laboratory	24	Viscometers	45
Homogenisers	48	Vortexers	44
Hotplates	33	Water Baths - Unstirred	46
Ice Makers	15	Water Purification	22
Incubators	8		

Over 300 special offers inside



Ultra-low Temperature Freezers

Haier Biomedical UK

Up to **46%** Off



UltraECO ULT Freezers Ranked #1 in ENERGY STAR RANKING



Energy-efficient design

- The UltraECO has the lowest energy consumption in the ENERGY STAR ranking – just 4.16 kWh/day. Ultra-low energy for ultra-low impact, making it a sustainable choice.

Save up to thousands of pounds

- You could save thousands of pounds in energy costs over 10 years with the UltraECO. We can also calculate your exact savings compared to your current unit.

Safety and reliability

- Dual sensors and self-diagnostic features ensure precise temperature control and quick alerts to prevent sample loss. VIP insulation and optional backup systems protect samples during both short and extended power outages.

IN STOCK NOW



Code	Type	Capacity (l)	Energy consumption (kwh/day)	List Price	Offer Price
DW-86L580BPT	Upright	580	5.18	£16,926.00	£9,795.00
DW-86L730BPT	Upright	730	4.16	£18,642.00	£9,995.00
DW-86L830BPT	Upright	830	5.18	£19,553.00	£10,678.00

J Series ULT Freezers

The Haier Biomedical J Series freezers deliver reliable ultra-low temperature performance. Designed for durability and ease of use, ensure safe, long-term storage of valuable samples.



DW-86L51J

DW-86L100J

DW-86L338J

DW-86L579


Up to **43%** Off

Haier Biomedical UK

Code	Type	Capacity (l)	Energy consumption (kwh/day)	List Price	Offer Price
DW-86L51J	Upright	51	3.7	£6,426.00	£4,347.00
DW-86L100J	Upright	100	5.5	£9,075.00	£6,113.00
DW-86L338J	Upright	338	7.5	£11,179.00	£6,345.00
DW-86L419	Upright	419	8	£11,474.00	£6,995.00
DW-86L579	Upright	579	9	£13,711.00	£8,295.00

Binder UFV 500 Ultra-Low Temperature Freezer

Code	Type	Capacity (l)	Energy consumption (kwh/day)	List Price	Offer Price
9020-0347	Upright	475	7.9	£13,156.00	£6,995.00



Class II Cabinets

Haier Biomedical UK

Haier Biomedical Class II Cabinets

Haier Biomedical Class II cabinets provide reliable personnel, product, and environmental protection for everyday laboratory work. Engineered for safety and efficiency, they deliver consistent performance across demanding applications.

Key Features & Benefits:

- Advanced airflow system for stable protection and reduced contamination risk
- Low energy consumption to cut running costs and support sustainability
- Quiet operation for a more comfortable working environment
- Intuitive touch-screen control panel for simple monitoring and adjustments
- Robust construction ensuring long-term durability and reliability



Code	Front closure	Working width (mm)	Noise level (dB)	Options included	List Price	Offer Price
HR900-IIA2-D-NS	Manual sliding sash	920	60	UV lamp, double electric socket	£9,878.00	£6,200.00
HR1200-IIA2-D-NS	Manual sliding sash	1,310	58	UV lamp, double electric socket	£10,656.00	£6,561.00

Haier Biomedical UK

Up to **49%** Off

Haier Biomedical Bio-2T Portable Liquid Nitrogen Sample Storage

The Bio-2T Portable is a compact, lightweight LN₂ transfer solution designed for the safe, short-distance transport of valuable biological samples in laboratory and clinical environments.

Key Features & Benefits:

- Lightweight sample transfer vessel
- 54 x 2ml Vial storage
- Real time temperature display
- High thermal insulation performance
- Anti-spill design for enhanced safety



Code	LN2 Capacity	Type	Capacity of 2ml vials	Autofill facility	List Price	Offer Price
Bio-2T	2	Tube	54	No	£1,393.00	£710.00

Up to **51%** Off

Haier Biomedical YDC-300H CryoStation

The Haier Biomedical YDC-3000H CryoStation is a mobile liquid nitrogen platform designed for the safe, controlled transport and handling of biological samples at cryogenic temperatures down to -196°C.

Key Features & Benefits:

- Real-time monitoring system
- LCD touchscreen control
- USB data export allows secure temperature data recording to support compliance and audits.
- Smooth-rolling, durable build for safe and efficient transport around the lab.



Code	LN2 Capacity	Type	Capacity of 2ml vials	Autofill facility	List Price	Offer Price
YDC-3000H	33	Box	50	No	£18,395.00	£9,043.00

Centrifuges & Microfuges



Up to **45%** Off

Save More with our Hettich Centrifuge Packages

Get accessories tailored to your application while saving more compared to purchasing units and accessories separately. We offer a wide range of package deals across our website—here are some of our favourites, all available and IN STOCK now.



Each package includes a swing-out rotor, buckets, lids, and adapters for the relevant tubes.



IN STOCK NOW



Code	Model	Type	Max g force	Application	Refrigerated	List Price	Offer Price
GB-C1406-P1	Universal 320R	Benchtop	21,382	Tissue culture - low capacity	Yes	£9,639.00	£6,073.00
GB-C1206-P5	Rotofix 32A	Benchtop	4,186	Blood tubes - low capacity	No	£5,251.00	£3,435.00
GB-C1706-P1	Rotina 380R	Benchtop	24,400	Tissue culture - medium capacity	Yes	£13,121.00	£8,582.00
GB-C1203-P1	Mikro 185	Microfuge	17,008	24 x 2ml tubes	No	£2,346.00	£1,431.00
GB-C2405-P1	Mikro 200R	Microfuge	21,382	24 x 2ml tubes	Yes	£6,504.00	£4,257.00

Shop our Hettich Centrifuges also IN STOCK

Code	Model	Type	Max g force	Refrigerated	List Price	Offer Price
GB-C1206	Rotofix 32A	Benchtop	4,186	No	£2,725.00	£1,795.00
GB-C4701	Rotofix 32A	Benchtop	24,400	No	£6,543.00	£3,615.00
GB-C1406	Universal 320R	Benchtop	21,382	Yes	£7,891.00	£4,609.00
GB-C5660	Rotanta 460R	Benchtop	24,400	Yes	£11,665.00	£7,619.00

OHAUS Benchtop Centrifuge Packages

Streamline your workflow with complete FC benchtop centrifuge packages, combining high-performance models with perfectly matched rotors and accessories for efficient, ready-to-run cell culture and sample processing.



Code	Model	Type	Max g force	Refrigerated	Capacity	List Price	Offer Price
30314812-Package2	FC5718 +	Benchtop	23,542	No	8 x 15ml or 4 x 50ml	£4,682.00	£2,945.00
83042045-Package2	FC5718R +	Benchtop	23,542	Yes	8 x 15ml or 4 x 50ml	£7,611.00	£4,787.00
30314816-Package1	FC5816 +	Benchtop	21,379	No	28 x 15ml or 12 x 50ml	£6,227.00	£3,917.00
83042047-Package2	FC5816R +	Benchtop	24,325	Yes	28 x 15ml or 12 x 50ml	£8,976.00	£5,584.00
30553036-Package1	FC5916 +	Benchtop	24,325	No	68 x 15ml or 28 x 50ml	£8,966.00	£5,640.00
83042049-Package3	FC5916R +	Benchtop	26,331	Yes	68 x 15ml or 28 x 50ml	£11,905.00	£7,488.00
83042055-Package1	FC5917RF +	Benchtop	26,361	Yes	48 x 15ml or 20 x 50ml	£14,472.00	£9,102.00

DLab D1008 Palm Microfuge

The DLab D1008 Palm Microfuge is a compact, easy-to-use microcentrifuge designed for rapid spin-downs and routine sample processing in busy laboratory environments.



Code	Model	Max g force	List Price	Offer Price
9031002012	D1008 Palm Green	7,000	£168.00	£99.00

IN STOCK NOW

CappRondo Mini Microfuge



Code	Max g force	Model	List Price	Offer Price
CR-1512	15,596	CAPPRondo Mini	£1,383.00	£694.00

IN STOCK NOW

CO₂ Incubators



Best conditions for your success

Up to **45%** Off

Binder CB-170 Solid Line CO₂ Incubator

The CB170 Binder CO₂ incubator delivers precise temperature and gas control to ensure optimal cell culture conditions. Its reliable performance and uniform environment supports consistent, reproducible results.

Key Features & Benefits:

- Fail-Safe CO₂ System: Automatically switches to a backup CO₂ supply if the primary source fails, ensuring continuous CO₂ control whilst protecting valuable cell cultures.
- Stable CO₂ and temperature regulation
- Seamless stainless steel inner chamber
- Humidity system for consistent conditions

Fail-Safe CO₂ System



Code	Model	Volume (l)	Disinfection routine	CO ₂ System Fail-Safe System	List Price	Offer Price
9640-0001	CB-S170 Solid Line	170	Heat 180°C	No	£7,625.00	£4,160.00
9640-0009	CB170	170	Heat 180°C	Yes	£10,357.00	£5,601.00
9040-0207	CB260	260	Heat 180°C	Yes	£13,916.00	£9,275.00



IN STOCK NOW 

9640-0001 & 9640-0009 MODELS ARE IN STOCK

GS Biotech 170LGOLD/HT CO₂ Incubator

Code	Volume (l)	Disinfection routine	List Price	Offer Price
170LGOLD/HT	170	Heat 120°C	£6,458.00	£3,995.00

IN STOCK NOW 



Cooled Incubators



Up to **21%** Off

LMS Series 4 600NP Cooled Incubator

The Series 4 600NP cooled incubator delivers exceptional precision temperature control in a large-capacity 600L stainless steel cabinet, making it ideal for demanding research and production environments. Designed for flexibility and reliability, it supports complex incubation protocols while maintaining outstanding temperature stability across the chamber.

Key features and benefits:

- **Fully programmable control (NP model)** – Enables multi-step temperature and time cycling for advanced experimental protocols.
- **Wide temperature range** – Operates from -10°C to +50°C for versatile application use.
- **Fan-assisted air circulation** – Ensures consistent and uniform temperature distribution throughout the chamber.
- **Robust stainless steel construction** – Provides durability, easy cleaning, and long-term performance in laboratory environments.



Code	Capacity (l)	Programmable controller	List Price	Offer Price
600NP	600	Yes	£10,812.00	£8,592.00

Incubators



Up to **27%** Off

Best conditions for your success

Binder BF115 Avantgarde.Line and BD-S56 Solid.Line Incubators

The Binder BD-S56 Solid.Line incubator is a simple, cost-effective unit designed for routine applications, offering reliable temperature control with natural convection and a straightforward interface. In contrast, the Binder BF115 Avantgarde.Line provides enhanced performance with forced-air circulation for improved uniformity, greater capacity, and advanced features such as high-temperature disinfection, making it better suited to demanding, precision-driven laboratory work.



Code	Model	Volume (l)	Fan	List Price	Offer Price
9010-0315	BF115 Avantgarde.Line	114	Yes	£2,769.00	£2,011.00
9090-0016	BD-S56 Solid.Line	56	Yes	£1,223.00	£889.00

Shaking Incubators



Up to **39%** Off

Be Scientific INS Shaking Incubators

The Be Scientific benchtop shaking incubators combine precise temperature control with reliable orbital shaking performance, making them ideal for applications including cell culture, microbiology and biochemical research. Available in a choice of capacities, these versatile units deliver smooth, uniform shaking alongside intuitive digital controls.

Key Features & Benefits:

- Adjustable Shaking Speed 25-600rpm
- Temperature range 5°C above room temperature to 65°C on INS-525 & INS-4100 models. INS-4200R can operate between 4-70°C
- Microprocessor control with large LCD screen, showing set and actual parameters (INS-525 & INS-4100)



Code	Refrigerated	Max vessel capacity (ml)	Max shaking capacity (ml)	List Price	Offer Price
INS-525	No	250	1,500	£2,921.00	£1,801.00
INS-4100	No	1,000	5,000	£3,778.00	£2,301.00
INS-4200R	Yes	2,000	8,000	£5,887.00	£3,617.00

Autoclaves



Up to **28%** Off

Astell Benchtop and Compact Autoclaves

The Astell AMA compact autoclaves feature a top-loading vertical chamber design ideal for taller items and laboratories where floor space is limited, while the AMB benchtop range offers a front-loading configuration better suited to bench installation and easier loading of trays and instruments.

Key Features & Benefits:

- Fully programmable 5.7" colour touchscreen controller with USB data port
- Plug & Play Installation
- Ecofill models feature an integrated reservoir and water recirculation system, reducing manual filling and supporting up to 15 cycles between refills.



Packages Include Load Sensed Process Timing and Fan Cooling (External)

Code	Model	Capacity (l)	Loading	Autofill	List Price	Offer Price
AMB230-SC	AMB230 Benchtop Ecofill package	43	Front	Yes	£12,784.00	£9,553.00
AMB420-SC	AMB420 Classic package	33	Front	No	£10,274.00	£7,318.00
AMB430-SC	AMB430 Classic package	43	Front	No	£11,695.00	£8,485.00
AMA240-SC	AMA240 Compact Ecofill package	63	Top	Yes	£14,750.00	£11,101.00

Rodwell autoclaves are UK-manufactured using premium materials and advanced technology, offering custom-built solutions with reliable performance, low running costs, and intelligent diagnostics, all supported by a nationwide service network.

Rodwell Portaclave Autoclaves

The Ambassador is Rodwell's most popular mid-sized autoclave, available in 100L and 158L formats. It offers flexible configurations for all load types, with in-chamber heating or integrated steam generation, plus optional pre- and post-vacuum cycles and a range of loading accessories.

Key Features & Benefits:

- UK-manufactured with a strong focus on quality and longevity
- Energy-efficient design to help reduce running costs
- Quiet operation, ideal for busy lab environments
- Optional array of furniture such as loading baskets, waste containers, loading trollies and transfer carriages.



Code	Model	Capacity (l)	Loading	Autofill	List Price	Offer Price
ST1528E	Portaclave 15E	15	Top	No	£3,233.00	£2,458.00
ST2228E	Portaclave 22E	22	Top	No	£3,864.00	£2,926.00
ST1528V	Portaclave 15V	15	Top	No	£3,965.00	£3,000.00
ST3028E	Portaclave 30E	30	Top	No	£4,469.00	£3,375.00
ST2228V	Portaclave 22V	22	Top	No	£4,595.00	£3,468.00
ST3028V	Portaclave 30V	30	Top	No	£5,220.00	£3,933.00

Prestige Medical Autoclaves

Code	Model	Capacity (l)	Loading	Autofill	List Price	Offer Price
210048	Classic Media Extendable	12	Top	No	£1,485.00	£995.00
210001	Classic 126	9	Top	No	£1,269.00	£845.00
210004	Classic Gauged	12	Top	No	£1,602.00	£1,095.00

IN STOCK NOW



Excelsior Scientific

Excelsior Scientific Autoclave Indicator Test Strips

Code	Model	Sterilisation criteria	List Price	Offer Price
TST-4SE	Good to Go TST-4SE	Steam steriliser operating 121°C to 134°C	£68.00	£55.00
TST-6SE	Good to Go TST-6SE	Steam steriliser operating 121°C to 134°C	£73.00	£66.00
TST-5SE	Good to Go TST-5SE	Steam steriliser operating 121°C to 134°C	£73.00	£66.00



Freezers -20°C

Haier Biomedical UK

Up to **37%** Off

Haier Biomedical DW-25L92SF -20°C Freezer

The Haier DW-25L92SF -20°C Freezer is a compact and dependable laboratory freezer designed to provide stable and consistent temperature control for the secure storage of temperature-sensitive samples.

Key Features & Benefits:

- Compact -20°C laboratory freezer
- Reliable and consistent temperature stability
- Ideal for storing temperature-sensitive samples
- Space-saving design for busy labs
- Easy-to-use digital temperature control



Code	Capacity (l)	Safety level	Lockable	List Price	Offer Price
DW-25L92SF	92	ATEX	Yes	£1,476.00	£931.00



Up to **30%** Off

LEC LSFSF110 and LSFSF312BT-UK -20°C Freezers

The LEC LSFSF110 and LSFSF312BT-UK -20°C freezers continue to be two of LEC's most popular laboratory freezers, offering laboratories a reliable and cost-effective solution without compromising on performance. Designed to deliver dependable temperature stability, practical storage and low running costs.

Key Features & Benefits:

- Spark-free interior design to reduce ignition risk when storing sensitive materials
- Energy-efficient operation helping to lower day-to-day running costs
- Lockable doors supplied with keys for added security and controlled access
- Reliable temperature performance with audible and visual alarm protection



Code	Model	Capacity (l)	Safety level	Lockable	List Price	Offer Price
444411345	LSFSF110 Laboratory Essential	98	Sparkfree	Yes	£412.00	£290.00
444411366	LSFSF312BT-UK Laboratory Plus	312	Sprakfree	Yes	£1,543.00	£1,082.00

Haier Biomedical UK

Up to **37%** Off

Haier Biomedical HLR-310SF Fridge

The Haier HLR-310SF Laboratory Fridge is a reliable and efficient refrigeration solution designed to maintain stable and uniform temperatures for the safe storage of temperature-sensitive samples and reagents.

Key Features & Benefits:

- Stable and uniform temperature control
- Reliable storage for samples and reagents
- Energy-efficient operation
- User-friendly digital controller
- Compact footprint with generous capacity



Code	Capacity (l)	Safety level	Lockable	List Price	Offer Price
HLR-118SF	118	ATEX certified	Yes	£1,476.00	£931.00
HLR-310SF	310	ATEX certified	Yes	£1,878.00	£1,177.00



Up to **30%** Off

LEC LSFSR110 and LSFSR310BT Fridges

The LSFSR310BT and LSFSR110 are two of LEC's most popular refrigeration solutions, offering reliable laboratory-grade performance at a highly cost-efficient price point. Designed to combine dependable temperature control with practical everyday usability, both units provide excellent value for laboratories looking for trusted storage without the premium cost.

Key Features & Benefits:

- Spark-free interior for safe storage of sensitive materials
- Flexible shelving positions for adaptable storage
- Lockable door supplied with two keys
- Energy-efficient performance that delivers reduced running costs over time



Code	Model	Capacity (l)	Safety level	Lockable	List Price	Offer Price
444411344	LSFSR110 Laboratory Essential	158	Sparkfree	Yes	£412.00	£290.00
444411360	LSFSR310BT Laboratory Plus	310	Sparkfree	Yes	£1,357.00	£946.00

Fridge Freezers



Up to **30%** Off



LEC LSFSC138DC-UK and LSFSC280BT Fridge Freezers

Designed for reliable laboratory cold storage, the LSFSC138DC-UK and LSFSC280BT combine secure sample protection with practical everyday performance. Featuring spark-free interiors, lockable doors, access ports for monitoring probes, and flexible storage configurations, these units are ideal for laboratories requiring dependable temperature-controlled storage. With compact undercounter and larger upright options available, the range offers a solution to suit a variety of laboratory environments.



Code	Model	Fridge Capacity (l)	Freezer Capacity (l)	List Price	Offer Price
444411368	LSFSC138DC-UK Laboratory Essential	118	18	£853.00	£592.00
444411369	LSFSC280BT Laboratory Plus	182	98	£1,388.00	£1,005.00



Up to **30%** Off



Labcold RLPR0517-RLVF0417-STACKINGKITDUAL Fridge Freezer

The Labcold RLPR0517 fridge, RLVF0417 freezer and STACKINGKITDUAL offer a flexible, space-saving refrigeration solution that is extremely popular in laboratory environments. Combining reliable performance with a highly cost-efficient way to achieve dual-temperature storage without the expense of a single integrated unit.

Key Features & Benefits:

- Space-saving design combines fridge and freezer vertically to maximise floor space
- Fridge and freezer run separately for improved reliability and easier maintenance
- Delivers dual-temperature storage at a lower purchase and lifecycle cost than traditional fridge/freezers

Code	Fridge Capacity (l)	Freezer Capacity (l)	Safety level	Lockable	List Price	Offer Price
RLVF0417-RLPR0517-Package1	150	124	Sparkfree	Yes	£2,509.00	£1,746.00

Pharmacy Fridges



Up to **24%** Off

Labcold Pharmacy Fridges

The Labcold IntelliCold® Pharmacy Fridge range delivers advanced, compliant temperature-controlled storage designed specifically for the safe preservation of vaccines and temperature-sensitive medicines in clinical and pharmacy environments.

Key Features & Benefits:

- IntelliCold® touch-screen controller with intuitive operation and secure access
- 24/7 automated temperature logging with Bluetooth connectivity and data access
- Advanced alarm system including high/low temperature, door open and power failure
- 5 year warranty



Code	Capacity (l)	Lock	Solid/glass door	List Price	Offer Price
RLDF0524	150	Key	Solid	£1,427.00	£1,083.00
RLDG0524	150	Key	Glass	£1,790.00	£1,358.00
RLDF1024	340	Key	Solid	£1,908.00	£1,447.00
RLDG1024	340	Key	Glass	£2,489.00	£1,888.00
RLDF1524	440	Key	Solid	£3,377.00	£2,561.00
RLDG1524	440	Key	Glass	£3,634.00	£2,756.00



Up to **25%** Off

Dulas DLS145PG Pharmacy Fridge

Code	Capacity (l)	Lock	Solid/glass door	List Price	Offer Price
DLS145PG	145	Key	Glass	£1,891.00	£1,423.00



Scotsman®

Up to **49%** Off

Scotsman Ice Makers

Reliable, efficient, and built to last — Scotsman ice makers deliver consistent, high-quality ice for any commercial setting. With robust construction and intuitive controls, they provide dependable performance day after day.

Key Features & Benefits:

- Produces crystal-clear, consistent ice
- Energy-efficient design for lower running costs
- Durable stainless-steel build for long service life
- Intelligent control system for easy operation
- Wide range of capacities to suit any workspace



Code	Output (per day) (kg)	Bin capacity (kg)	Type	Type of ice	List Price	Offer Price
AF87AX-6	70	25	Underbench	Flake	£5,776.00	£2,950.00
AF107AX-6	130	30	Free standing	Flake	£7,182.00	£3,910.00
EF107AX-6	130	30	Free standing	Flake	£7,728.00	£4,187.00
AF127AX-6	126	40	Underbench	Flake	£7,784.00	£4,227.00
EF127AX-6	126	40	Underbench	Flake	£8,293.00	£4,183.00
AF157AX-6	150	60	Free standing	Flake	£9,352.00	£5,019.00
EF157AX-6	150	60	Free standing	Flake	£9,890.00	£5,292.00



EF127AX-6



AF107AX-6



All Scotsman Ice Makers feature XSafe, a natural sanitation system that runs 24/7, killing over 99% of viruses and bacteria (including coronaviruses) and helping prevent mould, slime, and odours.

Service and Calibration

Airflow Cabinets (Class I/Class II/Class III)

Autoclaves

Balances

Centrifuges

Climate Stability Chambers

Cryogenic Storage

Freeze Dryers

Fridges/Freezers/Cold Rooms

Ice Makers

Incubators

Osmometers

Pipettors

Temperature Monitoring Systems

Vacuum Ovens

Water Baths

Water Chillers

Water Purification Systems

... and many more

We offer UK wide service and calibration on leading brand names including:



Enhanced Service Support Through Our Partnership with



To strengthen our service offering, we have partnered with Labserve, one of the UK's leading independent laboratory equipment service providers. This gives our customers access to a nationwide network of specialist engineers, ISO 17025 accredited calibration, maintenance and repair services, and expert support across a wide range of laboratory equipment. Together, we deliver responsive, compliant and reliable solutions that maximise equipment performance and minimise downtime.

KEEP YOUR LAB RUNNING WITHOUT INTERRUPTION

Whether you need urgent repairs, routine servicing, calibration, or validation, our experienced engineers are ready to help.

Speak to our Service Team

✉ service@wolflabs.co.uk

☎ 01759 301142 (Option 3)



Analytical Balances



Up to **32%** Off

Kern & Sohn ABP Analytical Balances

Premium KERN ABP analytical balances deliver ultra-fast, highly stable weighing performance with advanced Single-Cell technology for exceptional precision in demanding laboratory applications. Featuring intuitive operation, automatic internal adjustment and high-quality draught shield protection, the ABP range is designed for reliable, efficient daily use.



Code	Max capacity (g)	Readability (g)	Type of Calibration	List Price	Offer Price
ABP100-4M	120	0.0001	Internal	£1,956.00	£1,328.00
ABP200-4M	220	0.0001	Internal	£2,401.00	£1,630.00
ABP300-4M	220	0.0001	Internal	£2,667.00	£1,811.00
ABP100-5M	135	0.00001	Internal	£3,734.00	£2,535.00
ABP200-4AM	220	0.0001	Internal	£5,156.00	£3,500.00
ABP200-5M	220	0.00001	Internal	£5,334.00	£3,621.00
ABP300-4AM	320	0.0001	Internal	£5,601.00	£3,802.00
ABP100-5DAM	120	0.0001	Internal	£5,867.00	£3,983.00
ABP200-5DAM	220	0.0001	Internal	£7,823.00	£5,310.00
ABP100-5AM	135	0.00001	Internal	£8,001.00	£5,431.00
ABP100-5DM	120	0.0001	Internal	£2,845.00	£1,932.00
ABP200-5DM	220	0.0001	Internal	£3,023.00	£2,052.00
ABP200-5AM	220	0.00001	Internal	£9,778.00	£6,638.00

Up to **32%** Off

Kern & Sohn YKT-01 Needle Printer

Complete your analytical balance setup with the KERN YKT-01 printer – ideal for creating clear, GLP-compliant hard copy records of weighing results, including date, time and sample information. With reliable dot matrix printing and advanced documentation functions, it's perfect for traceable laboratory workflows and quality-controlled environments



Code	List Price	Offer Price
YKT-01	£516.00	£351.00

Kern & Sohn ADT Analytical Balances

The Kern & Sohn ADT analytical balances combine exceptional weighing precision with automatic internal calibration and advanced data connectivity, ensuring reliable, compliant results while simplifying daily laboratory workflows

Key Features & Benefits:

- **Automatic Internal Calibration** – Maintains weighing accuracy automatically, reducing user intervention and ensuring consistent, reliable results.
- **Modern Touch Display** – Intuitive touchscreen operation enables quick access to functions, improving efficiency and ease of use.
- **GLP/ISO & FDA 21 CFR Part 11 Support** – Enhances traceability and data integrity, helping laboratories meet regulatory and quality requirements.
- **USB & RS-232 Connectivity with Windows Direct** – Transfers weighing data directly to PCs and laboratory systems, reducing manual data entry and transcription errors.



Code	Max capacity (g)	Readability (g)	List Price	Offer Price
ADT100-4	120	0.0001	£1,512.00	£1,027.00
ADT200-4	220	0.0001	£1,689.00	£1,148.00
ADT300-4	320	0.0001	£1,925.00	£1,315.00

Kern & Sohn ADS Analytical Balances

The Kern & Sohn ADS analytical balances deliver exceptional weighing precision and dependable performance, making them an ideal choice for routine analytical, research and quality control applications. Featuring external adjustment, the ADS range allows calibration using certified external weights, giving laboratories greater flexibility, full traceability and enhanced control over calibration procedures while maintaining a cost-effective solution for regulated environments.



Code	Max capacity (g)	Readability (g)	List Price	Offer Price
ADS100-4	120	0.0001	£1,436.00	£984.00
ADS200-4	220	0.0001	£1,614.00	£1,104.00
ADS300-4	320	0.0001	£1,791.00	£1,225.00

Sartorius Entris II Essential BCE Analytical Balances

Delivering reliable precision and everyday ease of use, the Sartorius Entris® II Essential Line balances combine advanced weighing technology, intuitive operation and dependable performance for routine laboratory applications.

Key Features & Benefits:

- High-accuracy weighing for reliable laboratory results
- Fast stabilisation times to improve workflow efficiency
- Easy-to-use touchscreen display for simple operation
- Compact, durable design ideal for everyday lab use



Code	Max capacity (g)	Readability (g)	Type of Calibration	List Price	Offer Price
BCE64-1S	60	0.0001	External	£1,830.00	£1,057.00
BCE64I-1S	60	0.0001	Internal	£2,031.00	£1,172.00
BCE124-1S	124	0.0001	External	£1,951.00	£1,126.00
BCE124I-1S	120	0.0001	Internal	£2,142.00	£1,236.00
BCE224-1S	220	0.0001	External	£2,082.00	£1,201.00
BCE224I-1S	220	0.0001	Internal	£2,243.00	£1,249.00

Microbalances



Kern & Sohn ABT Microbalances

The Kern & Sohn ABT Microbalance range delivers exceptional weighing accuracy down to 0.01 mg, combining advanced single-cell technology, automatic internal calibration and rapid stabilisation times for reliable, high-precision laboratory measurements.



Code	Max capacity (g)	Readability (mg)	Dual or triple capacity	List Price	Offer Price
ABT100-5NM	101	0.01	No	£2,178.00	£1,479.00
ABT220-5DNM	220	0.1	Dual	£2,312.00	£1,570.00

Sartorius Entris II Essential BCE Top Pan Balances

Reliable, durable and designed for long-term performance, the Sartorius BCE Top Pan Balances deliver precise, repeatable results in both teaching environments and professional laboratories. With intuitive operation, a clear backlit display and a tough, chemical-resistant housing, they're built to stand up to daily use – from routine weighing to more demanding applications.

Key Features & Benefits:

- **High accuracy and broad capacity range** – ideal for diverse sample types and tasks.
- **Fast stabilisation** – streamline workflows and minimise waiting time between measurements.
- **Robust, easy-clean construction with overload protection** – dependable in classroom, research or industrial settings.
- **Flexible connectivity (USB, RS232)** – supports efficient data transfer and integration with external systems.



Code	Max capacity (g)	Readability (g)	Type of calibration	List Price	Offer Price
BCE223-1S	220	0.001	External	£1,135.00	£661.00
BCE223I-1S	220	0.001	Internal	£1,367.00	£793.00
BCE323-1S	320	0.001	External	£1,246.00	£724.00
BCE323I-1S	320	0.001	Internal	£1,477.00	£856.00
BCE423-1S	420	0.001	External	£1,377.00	£799.00
BCE423I-1S	420	0.001	Internal	£1,578.00	£914.00
BCE622-1S	620	0.01	External	£1,027.00	£599.00
BCE622I-1S	620	0.01	Internal	£1,236.00	£719.00
BCE623-1S	620	0.001	External	£1,608.00	£931.00
BCE623I-1S	620	0.001	Internal	£1,800.00	£1,040.00
BCE653-1S	650	0.001	External	£1,830.00	£1,057.00
BCE653I-1S	650	0.001	Internal	£2,031.00	£1,172.00
BCE822-1S	820	0.01	External	£1,135.00	£661.00
BCE822I-1S	820	0.01	Internal	£1,367.00	£793.00
BCE1202-1S	1,200	0.01	External	£1,246.00	£724.00
BCE1202I-1S	1,200	0.01	Internal	£1,477.00	£856.00
BCE2201-1S	2,200	0.1	External	£1,135.00	£661.00
BCE2201I-1S	2,200	0.1	Internal	£1,367.00	£793.00
BCE2202-1S	2,200	0.01	External	£1,377.00	£799.00
BCE2202I-1S	2,200	0.01	Internal	£1,578.00	£914.00
BCE3202-1S	3,200	0.01	External	£1,488.00	£862.00
BCE3202I-1S	3,200	0.01	Internal	£1,709.00	£989.00
BCE4202-1S	4,200	0.01	External	£1,608.00	£931.00
BCE4202I-1S	4,200	0.01	Internal	£1,800.00	£1,040.00
BCE5201-1S	5,200	0.1	External	£1,246.00	£724.00
BCE5201I-1S	5,200	0.1	Internal	£1,477.00	£856.00
BCE6200-1S	6,200	0.01	External	£1,135.00	£661.00
BCE6200I-1S	6,200	1	Internal	£1,367.00	£793.00
BCE6202-1S	6,200	0.1	External	£1,830.00	£1,057.00
BCE6202I-1S	6,200	0.01	Internal	£2,031.00	£1,172.00
BCE8200-1S	8,200	1	External	£1,246.00	£724.00
BCE8200I-1S	8,200	1	Internal	£1,477.00	£856.00
BCE8201-1S	8,200	0.1	External	£1,377.00	£799.00
BCE8201I-1S	8,200	0.1	Internal	£1,578.00	£914.00

Portable Balances



Up to **35%** Off

Ohaus Navigator Portable Balances

The newly designed OHAUS Navigator series delivers excellent value for general weighing applications in laboratory, industrial and education environments. Rugged, multifunctional and easy to use, it features fast stabilisation, exceptional overload protection and reliable everyday performance.

Key Features & Benefits:

- Fast one-second stabilisation for rapid and efficient weighing results
- Superior overload protection capable of handling up to 400% of rated capacity
- Multiple weighing applications including parts counting, checkweighing and percent weighing
- Easy-to-read backlit LCD display for clear visibility in all working conditions



Code	Model	Max capacity (g)	Readability (g)	Type of calibration	List Price	Offer Price
30456385	Navigator NV221	220	0.1	External	£211.00	£140.00
30456392	Navigator NVT2200	2,200	1	External	£212.00	£141.00
30456393	Navigator NVT6200	6,200	1	External	£217.00	£144.00
30456386	Navigator NV621	620	0.1	External	£253.00	£166.00
30456387	Navigator NV1201	1,200	0.1	External	£285.00	£186.00
30456388	Navigator NV2201	2,200	0.1	External	£298.00	£195.00
30456389	Navigator NVT2201	2,200	0.1	External	£315.00	£205.00
30456394	Navigator NVT12000	12,000	1	External	£316.00	£206.00
30456390	Navigator NVT4201	4,200	0.1	External	£326.00	£211.00
30615581	Navigator NV123	120	0.001	External	£334.00	£217.00
30456382	Navigator NV222	220	0.01	External	£338.00	£220.00
30456391	Navigator NVT6201	6,200	0.1	External	£340.00	£221.00
30456383	Navigator NV422	420	0.01	External	£345.00	£224.00
30615587	Navigator NVT10201	10,200	0.1	External	£352.00	£228.00
30456384	Navigator NV622	620	0.01	External	£352.00	£228.00
30456395	Navigator NVT22000	22,000	1	External	£358.00	£232.00
30615584	Navigator NV1202	1,200	0.01	External	£358.00	£232.00
30615582	Navigator NV223	220	0.001	External	£370.00	£240.00
30615585	Navigator NV2202	2,200	0.01	External	£374.00	£242.00
30615583	Navigator NV323	320	0.001	External	£385.00	£349.00
30615586	Navigator NV3202	3,200	0.01	External	£394.00	£255.00

Water Purification Systems

WolfLabs works closely with Merck to provide high-quality water purification systems designed for reliable and consistent laboratory performance. Merck's systems deliver pure and ultrapure water solutions to meet the demands of a wide range of scientific applications.

Up to **47%** Off



Milli-Q SQ 2 Water Purification Packages

Water Quality

- Delivers ultrapure and pure water suitable for sensitive lab applications.
- Continuously monitors water quality to ensure consistent performance.

Design & Usability

- Compact, modern design saves valuable bench space.
- Intuitive touchscreen interface for easy operation and system monitoring.

Customisation & Efficiency

- Customisable water output tailored to specific lab needs.
- Intelligent system tracks water usage and predicts maintenance requirements.

Maintenance & Reliability

- Long-lasting purification cartridges reduce downtime.
- Low-maintenance design ensures consistent and reliable operation.



Code	Max output of type 1 water (l/min)	Max output of type 3 water (l/hr)	Plug and play?	List Price	Offer Price
ZSQ200UPT0-Package	1.6	0	Yes	£6,292.00	£3,329.00
ZSR240R0T0-Package	1.6	42	Yes	£7,319.00	£3,962.00
ZSQ240UCT0-Package	1.6	42	Yes	£7,581.00	£4,006.00



ZSQ200UPT0-Package
Pressurised loop configuration (suitable for labs with a pressurised purified water loop)



ZSR240R0T0-Package
Standard configuration (suitable for one small team)



ZSQ240UCT0-Package
Combined configuration (suitable for one small lab or team)

Recirculating Fume Cupboards



Up to **33%** Off

Esco Ascent Max Recirculating Fume Cupboards

Keep your team safe with clean, filtered air – without the hassle of ducting. The Ascent Max delivers high-performance fume containment in a plug-and-play design, saving space, time, and installation cost.

Key Features & Benefits:

- **Ductless design** – Easy setup, no building works
- **Advanced carbon filtration** – Helps protect users from harmful fumes
- **Sentinel™ control system** – Simple, safe, real-time monitoring
- **Large working area** – Comfortable, practical workspace
- **Front-access filter replacement** – Quick, minimal disruption
- **Durable steel construction** – Long-lasting reliability



We also have the required filters on special offer, contact us for a full quotation.

Code	Model	Nominal width (mm)	Low airflow alarm	Work surface included	List Price	Offer Price
2040018	Ascent Max ADC-2B1	600	Yes	Yes	£6,652.00	£4,479.00
2040024	Ascent Max ADC-3B1	900	Yes	Yes	£6,652.00	£4,479.00
2040033	Ascent Max ADC-3C1	1,035	Yes	Yes	£6,959.00	£4,678.00
2040039	Ascent Max ADC-3E1	900	Yes	Yes	£7,954.00	£5,323.00
2040042	Ascent Max ADC-4B1	1,200	Yes	Yes	£6,652.00	£4,479.00
2040051	Ascent Max ADC-4C1	1,340	Yes	Yes	£6,959.00	£4,678.00
2040057	Ascent Max ADC-4E1	1,200	Yes	Yes	£7,954.00	£5,323.00
2040060	Ascent Max ADC-5B1	1,500	Yes	Yes	£7,214.00	£4,843.00
2040069	Ascent Max ADC-5C1	1,645	Yes	Yes	£7,551.00	£5,061.00
2040075	Ascent Max ADC-5E1	1,500	Yes	Yes	£8,844.00	£5,899.00
2040078	Ascent Max ADC-6B1	1,800	Yes	Yes	£7,214.00	£4,843.00
2040087	Ascent Max ADC-6C1	1,950	Yes	Yes	£7,551.00	£5,061.00
2040213	Ascent Max ADC-6E1	1,800	Yes	Yes	£8,316.00	£5,448.00

Support Stands for ESCO Ascent Fume Cupboards

Code	Description	Compatible Fume Cupboards	List Price	Offer Price
5130169	Support Stand with Leveling Feet (Height 28") for 2ft cabinets	2040018	£766.00	£544.00
5130173	Support Stand with Leveling Feet (Height 34") for 2ft cabinets	2040018	£766.00	£544.00
5130170	Support Stand with Leveling Feet (Height 28") for 3ft cabinets	2040024 2040033 2040039	£766.00	£544.00
5130174	Support Stand with Leveling Feet (Height 34") for 3ft cabinets	2040024 2040033 2040039	£766.00	£544.00
5130134	Support Stand with Leveling Feet (Height 28") for 4ft cabinets	2040042 2040051 2040057	£766.00	£544.00
5130175	Support Stand with Leveling Feet (Height 34") for 4ft cabinets	2040042 2040051 2040057	£766.00	£544.00
5130171	Support Stand with Leveling Feet (Height 28") for 5ft cabinets	2040060 2040069 2040075	£817.00	£579.00
5130176	Support Stand with Leveling Feet (Height 34") for 5ft cabinets	2040060 2040069 2040075	£817.00	£579.00
5130172	Support Stand with Leveling Feet (Height 28") for 6ft cabinets	2040078 2040087 2040213	£817.00	£579.00
5130177	Support Stand with Leveling Feet (Height 34") for 6ft cabinets	2040078 2040087 2040213	£817.00	£579.00
5130161	Support Stand with Caster Wheels (Height 28") for 2ft cabinets	2040018	£680.00	£486.00
5130164	Support Stand with Caster Wheels (Height 34") for 2ft cabinets	2040018	£680.00	£486.00
5130155	Support Stand with Caster Wheels (Height 28") for 3ft cabinets	2040024 2040033 2040039	£680.00	£486.00
5130165	Support Stand with Caster Wheels (Height 34") for 3ft cabinets	2040024 2040033 2040039	£680.00	£486.00
5130152	Support Stand with Caster Wheels (Height 28") for 4ft cabinets	2040042 2040051 2040057	£680.00	£486.00
5130166	Support Stand with Caster Wheels (Height 34") for 4ft cabinets	2040042 2040051 2040057	£680.00	£486.00
5130162	Support Stand with Caster Wheels (Height 28") for 5ft cabinets	2040060 2040069 2040075	£764.00	£542.00
5130167	Support Stand with Caster Wheels (Height 34") for 5ft cabinets	2040060 2040069 2040075	£764.00	£542.00
5130154	Support Stand with Caster Wheels (Height 28") for 6ft cabinets	2040078 2040087 2040213	£764.00	£542.00
5130168	Support Stand with Caster Wheels (Height 34") for 6ft cabinets	2040078 2040087 2040213	£764.00	£542.00

Laboratory Glass Washers



Up to **35%** off

Miele ExploreLine Laboratory Glass Washers

The Miele ExploreLine Glasswasher is designed for reliable, everyday lab cleaning. It offers efficient wash programs, energy-saving operation, and consistent results, making it ideal for routine cleaning of standard labware.

Key Features & Benefits:

- **Efficient Wash Programs** – Thorough cleaning for standard labware with minimal water and energy use.
- **Consistent, Reproducible Results** – Reliable performance ensures labware is cleaned to the same high standard every time.
- **User-Friendly Controls** – Simple interface reduces training time and improves workflow efficiency.
- **Energy-Efficient Design** – Optimised cycles save power and water, lowering operational costs.



Packages include basic starter pack

Code	Model	Type	Max temp (°C)	Accessories included	List Price	Offer Price
12797660-Package1	ExploreLine PLW 8604 ADP	Under bench	70	2 baskets, 4 inserts	£8,979.00	£5,879.00
12797650-Package1	ExploreLine PLW 8604 AD	Under bench	70	2 baskets, 4 inserts	£8,705.00	£5,700.00



Miele ExpertLine Laboratory Glass Washers

The Miele ExpertLine Glasswasher is engineered for high-demand laboratories, offering advanced cleaning performance, faster cycles, and enhanced flexibility. It handles larger loads and a wider variety of labware while ensuring precise, reproducible results for even the most critical applications.

Key Features & Benefits:

- **Advanced Wash Programs** – Optimised for complex and heavy-duty labware cleaning.
- **High Capacity & Flexible Loading** – Accommodates larger volumes and diverse labware types.
- **Consistent, Reproducible Results** – Maintains precise cleaning standards for critical lab applications.
- **Fast Cycle Times** – Reduces downtime and improves lab productivity.

Packages include basic starter pack

Code	Model	Type	Max temp (°C)	Accessories included	List Price	Offer Price
12798760-Package1	ExpertLine PLW 8683 SST AD DC10	Under bench	93	2 baskets, 2 inserts, 2 modules	£14,439.00	£9,452.00




EXTRA

Upgrade Your Undercounter Glasswasher

For a limited time, upgrade and receive:

- FREE A101 upper basket & A150 lower basket
- Free WEEE collection & sustainable disposal of your old washer
- Connection to Miele MOVE Connect

Contact our Sales Team for details.



Up to **25%** Off

Introducing the Grant ORCA Recirculating Chiller

The ORCA Closed Loop Refrigerated Circulator delivers precise, energy-efficient cooling while eliminating water waste. With $\pm 0.1^\circ\text{C}$ stability and a compact, quiet design, it's a reliable, cost-effective upgrade for modern labs.

Temperature: -10 to 40°C operating range with superb $\pm 0.1^\circ\text{C}$ stability

New Tech: Advanced refrigeration automatically adjusts to heat loads, using only the power required

Reduce Waste: Replaces costly, water-intensive tap-water cooling systems

Ultra-Quiet Operation: Only 49 dB, ideal for lab environments.

Compact: Small footprint to minimise lost bench or floor space

Ease of Use: Intuitive, straightforward controls—no training needed

Modern Design: Showcases Grant's new, sleek brand identity - an attractive addition to any lab

Uncluttered Workspace: Rear-mounted connectors and switches keep the bench clear

Slashing Utilities: Dramatically lowers water use and recycles heat - reduction in water and HVAC costs



Code	Cooling capacity (kW)	Min temp ($^\circ\text{C}$)	Max temp ($^\circ\text{C}$)	Pump flow rate (l/min)	List Price	Offer Price
ORCA-UK	0.3	-10	40	15	£3,124.00	£2,333.00

Recirculating Chillers & Water Baths



Up to **33% off**

Grant ECOCOOL Recirculating Chiller

The Grant ECOCOOL recirculating chiller provides reliable, energy-efficient cooling for a wide variety of laboratory applications. Designed with performance, sustainability, and ease of use in mind, it delivers precise temperature control in a compact, quiet, and user-friendly system.

- **Energy-Efficient Cooling** – Designed to reduce running costs and environmental impact compared to traditional chillers.
- **Precise Temperature Control** – Delivers accurate and stable cooling for a wide range of laboratory applications.
- **Compact and Space-Saving** – Small footprint makes it ideal for crowded lab environments.
- **User-Friendly Interface** – Easy-to-use digital controls for quick setup and reliable operation.
- **Quiet Operation** – Low noise levels ensure minimal disruption in the laboratory.



Code	Cooling capacity (kW)	Min temp (°C)	Max temp (°C)	Pump flow rate (l/min)	List Price	Offer Price
ECOCOOL100R	0.25	-20	100	17	£4,112.00	£2,817.00
ECOCOOL150R	0.25	-30	150	17	£4,730.00	£3,332.00

Grant TXF200-R52 Recirculating Chiller

The Grant TXF200-R52 refrigerated circulating bath delivers ultra-precise temperature control from -52°C to 100°C, combining powerful 1.1kW cooling performance, ±0.1°C stability and programmable operation for demanding laboratory applications.



Code	Cooling capacity (kW)	Min temp (°C)	Max temp (°C)	Pump flow rate (l/min)	List Price	Offer Price
TXF200-R52	0.9	-52	100	23	£12,813.00	£8,899.00

Grant Unstirred Water Baths

Code	Tank capacity (l)	Control	Stability (+/- °C)	List Price	Offer Price
JBA5	5	Digital	0.5	£545.00	£371.00
JBN5	5	Digital	0.5	£698.00	£461.00
SAP2S	2	Digital	0.2	£783.00	£468.00
SAP5	2	Digital	0.2	£843.00	£552.00



IN STOCK NOW



Drying/Warming Ovens



Up to **10 %** Off

LEEC ECO Drying Cabinets

The LEEC ECO Drying Cabinets are designed to deliver reliable, energy-efficient drying performance for laboratory glassware, instruments and media. Manufactured in the UK, the range combines low power consumption with advanced thermal insulation to reduce heat loss and running costs, while maintaining excellent temperature stability.

Key Features & Benefits:

- Energy-saving thermal insulation for reduced heat loss and lower operating costs
- Digital controller with full colour touch screen display for easy operation
- Programmable 24/7 scheduling function to optimise energy usage
- Temperature range from 30°C to 80°C for versatile drying applications
- Anti-microbial powder coated exterior for added hygiene and durability



Code	Capacity (l)	Energy Efficient?	Overheat thermostat	List Price	Offer Price
ECOSS	113	Yes	Yes	£2,555.00	£2,288.00
ECOLS	227	Yes	Yes	£2,745.00	£2,456.00
ECOFC2	426	Yes	Yes	£4,265.00	£3,813.00
ECOFC1	885	Yes	Yes	£5,795.00	£5,176.00



Up to **20 %** Off

LTE ECO Drying Cabinets

The LTE Scientific ECO Drying Cabinets deliver an energy-efficient and sustainable solution for laboratory glassware drying, helping to reduce power consumption and running costs by up to 65% compared to traditional drying cabinets. Designed with five-sided insulation for improved heat retention, the range also features precise digital PID temperature control, a programmable 7-day timer and a durable stainless steel finish, combining reliable performance with everyday efficiency.



Code	Capacity (l)	Energy Efficient?	Overheat thermostat	List Price	Offer Price
TD0100SN	100	Yes	No	£2,049.00	£1,646.00
TD0180SN	180	Yes	No	£2,311.00	£1,848.00

Vacuum Ovens



Best conditions for your success

Up to **22%** Off

Binder VD and VDL Vacuum Ovens

The BINDER VD and VDL vacuum ovens are designed for fast, gentle and highly efficient drying of sensitive samples without residue, oxidation or contamination. The VD series is ideal for non-flammable solvents, while the ATEX-compliant VDL range provides maximum safety when working with flammable solvents.

Key Features & Benefits:

- **Gentle vacuum drying** – Protects sensitive samples and reduces drying times.
- **Uniform heating** – Excellent heat transfer for consistent results.
- **Enhanced safety** – ATEX-compliant VDL models for flammable solvents.
- **Easy operation** – Flexible shelving and intuitive controls.



Code	Model	Capacity (l)	Max temp (°C)	ATEX compliant	List Price	Offer Price
9630-0002	VD56	55	220	No	£8,444.00	£6,540.00
9630-0009	VDL23	24	110	Yes	£13,031.00	£10,094.00

General Purpose Vacuum Pumps



Up to **28%** Off

KNF LABOPORT SC Vacuum Pumps

The KNF LABOPORT® SC vacuum systems combine chemically resistant, oil-free performance with intuitive touchscreen vacuum control, delivering precise, quiet and environmentally friendly operation for demanding laboratory applications.



Code	Model	Type	Chemically resistant	Ultimate vacuum (mbar)	List Price	Offer Price
329849/338082	SC 820 G	Diaphragm	Yes	6	£5,426.00	£3,898.00
329851/338084	SC 840 G	Diaphragm	Yes	6	£5,915.00	£4,249.00
306473/117578	SC 920 G	Diaphragm	Yes	2	£5,184.00	£3,583.00

Syringe Pumps



Up to **25 %** off

Biochrom Legato Syringe Pumps

The KDS Legato Syringe Pumps are precision infusion and withdrawal devices designed for accurate, reliable fluid delivery in laboratory and research settings. Their compact, robust design makes them ideal for applications requiring consistent and reproducible flow rates.

Key Features & Benefits:

- **High-precision flow control** – delivers accurate and repeatable results.
- **Wide flow rate range** – supports both low and high flow applications.
- **Intuitive touchscreen interface** – enables quick setup and easy monitoring.
- **Broad syringe compatibility** – works with a wide range of syringe sizes.



Code	Max number of syringes	Max syringe volume (ml)	List Price	Offer Price
788100	1	60	£2,336.00	£1,735.00
788110	1	60	£3,045.00	£2,259.00

Spectrophotometers



Up to **39 %** off

DLab SP-UV1100 Spectrophotometer

The DLAB SP-UV1100 UV/Vis Spectrophotometer delivers reliable, high-precision analysis with a wide 190–1100nm wavelength range.

Key Features & Benefits:

- **190–1100nm UV/Vis wavelength range** – Suitable for a wide variety of analytical and research applications
- **2nm spectral bandwidth** – Provides accurate, reliable measurement results
- **Automatic lamp switching** – Seamlessly switches between tungsten and deuterium lamps for improved efficiency
- **Stores up to 200 data sets and curves** – Protects valuable data even after power interruption
- **Standalone or PC-controlled operation** – Flexible workflow with optional Wave Professional software



Code	Model	Min wavelength (nm)	Max wavelength (nm)	List Price	Offer Price
4010120102	SP-UV1100	190	1100	£3,943.00	£2,407.00

Overhead Stirrers



Up to **53 %** off

DLab OS20-S Overhead Stirrer Package

DLAB Overhead Stirrers are competent for most demanding tasks with guaranteed safety and performance. They are widely used in chemical synthesis, pharmaceutical synthesis, physical and chemical analysis, petrochemical industries, etc.

Key Features & Benefits:

- Wide range of speed from 50 to 2200rpm
- Brushless DC motor for long life, maintenance free and explosion proof
- LED screen for monitoring actual speed



Code	Max capacity	Speed range (rpm)	Package contents	List Price	Offer Price
8033100200+18900605	20	50-2200	Includes Universal plate stand (18900605) and Crossed stirrer (18900071)	£725.00	£337.00

IN STOCK NOW



Up to **27 %** off

IKA EUROSTAR 20 Digital Overhead Stirrer

The IKA EUROSTAR 20 digital overhead stirrer delivers precise, reliable mixing performance for a wide range of laboratory applications, combining quiet operation with consistent speed control for dependable results.

Key Features & Benefits:

- **Powerful stirring performance** – Handles viscosities up to 10,000 mPas with ease.
- **Wide speed range** – Adjustable speeds from 0/30 to 2,000 rpm for versatile applications.
- **Digital display** – Clear speed readout for accurate monitoring and control.
- **Quiet, compact design** – Space-saving operation ideal for everyday laboratory use.

Code	Max capacity	Speed range (rpm)	List Price	Offer Price
0020132602-IKA001	25	30-2000	£979.00	£715.00

Requires a stirrer attachment (sold separately)



Magnetic Stirrers

Cole-Parmer®

Up to **32% off**

Cole-Parmer ST-200 Magnetic Stirrers

Delivering reliable, high-performance mixing, the Cole-Parmer ST-200 Series stirrers combine powerful stirring capability with durable, space-saving designs for everyday laboratory applications.

Key Features & Benefits:

- **Powerful stirring performance** – Mix volumes up to 30 L with consistent, reliable speed control for demanding lab workflows.
- **Robust, durable construction** – Tough aluminium or stainless steel top plates provide long-lasting performance and chemical resistance.
- **Compact, space-saving design** – Ideal for busy laboratories with limited bench space and easy storage requirements.
- **Enhanced safety design** – Spill-resistant construction helps direct liquids away from controls and electronics.
- **BioCote® antimicrobial protection** – Built-in antimicrobial technology supports cleaner, more hygienic working environments



Code	Max capacity (l)	Control	Plate dimensions (w x d)	List Price	Offer Price
04807-54	15	Analogue	150 x 150mm	£351.00	£237.00
04807-74	30	Analogue	300 x 300mm	£948.00	£626.00

D LAB

Up to **36% off**

DLab EcoStir Magnetic Stirrers

EcoStir is a magnetic stirrer customized for liquid mixing, suitable for chemical analysis, biology reagent mixing, regular laboratory experiments and liquid handling.

Key Features & Benefits:

- Maintenance free brushless DC motor
- Rare earth magnets give powerful stirring
- Stepless speed regulation 300-2000rpm
- Ultra silence and vibration-free design
- PET top provides excellent chemical resistance
- Shaft can be fitted with a pH sensor or other sensors for different purposes



Code	Max capacity (l)	Control	Plate dimensions (dia)	List Price	Offer Price
8010194000	1.5	Analogue	120mm	£97.00	£62.00

IN STOCK NOW



IKA RH Basic, RH Lite and C-MAG HS 7 Magnetic Stirrers

Designed for reliable everyday laboratory mixing and heating applications, the IKA RH Basic, RH Lite and C-MAG HS 7 magnetic stirrers deliver precise temperature control, powerful stirring performance and durable construction for consistent, reproducible results.

Key Features & Benefits:

- Powerful stirring performance for reliable mixing of a wide range of samples
- Accurate temperature control for safe and consistent heating applications
- Durable chemical-resistant design ideal for demanding laboratory environments
- Easy-to-use controls with clear displays for simple day-to-day operation



Code	Max capacity (l)	Control	Plate dimensions (w x d or dia)	List Price	Offer Price
0010016456-IKA001	10	Digital	180 x 180mm	£969.00	£707.00
0025006148-IKA001	15	Analogue	135mm	£580.00	£426.00
0025006149-IKA001	15	Digital	135mm	£384.00	£285.00

C-MAG HS 7 Package includes **C-MAG HS 7** stirrer, **ETS-D5** Electronic contact thermometer, **H16 V** Support rod, **H38** Holding rod, and **H44** Boss head clamp

Hotplate Stirrers

DLab MS-H Hotplate Stirrers

Delivering reliable stirring performance with precise temperature control, the DLAB MS-H280-Pro and MS-H-ProA are ideal for routine laboratory mixing and heating applications, combining safety, efficiency and user-friendly operation in one compact solution.



8160221110 MS-H-ProA Package 1 includes PT1000A temperature sensor (**18900016**) and support clamp (**18900017**)

Code	Max capacity (l)	Max plate temp (°C)	Control	List Price	Offer Price
8030108212	3	280	Digital	£191.00	£119.00
8160221110	20	340	Digital	£586.00	£358.00

IN STOCK NOW



Hotplates

Cole-Parmer®

Up to **34%** Off

Cole-Parmer HP-200D Hotplates

The Cole-Parmer HP-200D large-capacity digital hot plates are designed to deliver reliable, precise heating for demanding laboratory applications. Featuring advanced microprocessor control and a spacious plate surface, these robust units are ideal for heating large vessels or multiple samples simultaneously while maintaining excellent temperature stability and user safety.

Key Features & Benefits:

- Accurate microprocessor-controlled heating for consistent temperature performance
- Large plate surface accommodates multiple vessels or oversized samples
- Hot warning indicator enhances user safety even after switch-off
- Built-in BioCote® antimicrobial protection for improved hygiene and cleaner operation



Code	Max plate temp (°C)	Control	Plate Material	List Price	Offer Price
04807-66	450	Analogue	Ceramic	£871.00	£596.00
04807-70	375	Analogue	Ceramic	£1,136.00	£755.00
04807-82	300	Digital	Aluminium	£1,194.00	£805.00

D LAB

Up to **39%** Off

D Lab HP-Pro Hotplate Packages

The DLAB HP-Pro Packages combine reliable heating performance with precise digital temperature control, making them ideal for a wide range of laboratory applications. Available in both the compact HP380-Pro and larger HP500-Pro models, these special offer packages include a PT1000A temperature sensor for enhanced accuracy and repeatable results. With durable ceramic-coated or glass ceramic plates, fast heat transfer, and intuitive digital controls, the HP-Pro range delivers dependable day-to-day performance while offering excellent value



Packages include temperature sensor **PT1000A**

Code	Max plate temp (°C)	Control	Plate Material	List Price	Offer Price
5131151114	380	Digital	Aluminium with ceramic coating	£307.00	£188.00
5131132115	500	Digital	Glass ceramic	£795.00	£485.00

IN STOCK NOW



Probe Type Sonicators



Up to **34%** Off

QSonica Probe Type Sonicators

The Q125, Q500 and Q700 Sonicators deliver powerful, reliable ultrasonic processing for sample preparation, cell disruption, homogenisation and nanoparticle applications.



Code	Max volume (l)	Wattage (W)	List Price	Offer Price
Q125AUK-220	0.05	125	£3,806.00	£2,501.00
Q500UK-220	1	500	£6,372.00	£4,200.00
Q700UK-220	1	700	£7,114.00	£4,674.00

Laboratory Chairs



Up to **26%** Off

PSI Seating PLAB Laboratory Chairs

Seating's PLAB laboratory chairs are designed specifically for scientific and technical environments, offering ergonomic support and easy-clean surfaces ideal for lab use. Their durable construction and adjustable features ensure comfort, safety, and long-lasting performance in demanding workspaces.



Castors can be swapped for glides – please notify our sales team if you'd like this option.

PLAB(1) PLAB(2) PLAB-STOOL PLAB-SMALL-STOOL

Code	Type	Material	Lift (height)	Glides or castors	List Price	Offer Price
PLAB(1)	Chair	Polyurethane	Bench	Castors	£165.00	£122.00
PLAB(2)	Chair	Polyurethane	Bench	Castors	£173.00	£128.00
PLAB-STOOL	Stool	Polyurethane	Desk	Castors	£147.00	£98.00
PLAB-SMALL-STOOL	Stool	Polyurethane	Desk	Castors	£134.00	£97.00

Manual Pipettes

IKA

Up to **56%** off

PETTE Vario



IKA PETTE Vario

The IKA PETTE vario offers precise, adjustable volume pipetting with excellent ergonomics for comfortable, reliable daily use. Its durable design and easy calibration make it ideal for accurate liquid handling across a wide range of applications.

PETTE Multi



IKA PETTE Multi

The IKA PETTE multi provides fast, efficient multichannel pipetting, perfect for plate-based workflows. With smooth operation and consistent dispensing across all channels, it's designed to maximise productivity while maintaining accuracy.

Code	Model	Max Number of channels	Volume range (min-max)	Autoclavable	List Price	Offer Price
0020011213-IKA001	PETTE Vario	1	2-20	Yes	£326.00	£144.00
0020011216-IKA001	PETTE Vario	1	100-1,000	Yes	£326.00	£144.00
0020011215-IKA001	PETTE Vario	1	20-2,000	Yes	£326.00	£144.00
0020011214-IKA001	PETTE Vario	1	10-100	Yes	£326.00	£144.00
0020011217-IKA001	PETTE Vario	1	500-5,000	Yes	£326.00	£144.00
0020011211-IKA001	PETTE Vario	1	0.5-10	Yes	£326.00	£144.00
0020011210-IKA001	PETTE Vario	1	0.1-2	Yes	£326.00	£144.00
0020011218-IKA001	PETTE Vario	1	1,000-10,000	Yes	£326.00	£144.00
0020017945-IKA001	PETTE Multi	8	20-200	Yes	£588.00	£266.00
0020017944-IKA001	PETTE Multi	8	10-100	Yes	£588.00	£266.00
0020017943-IKA001	PETTE Multi	8	0.5-10	Yes	£588.00	£266.00
0020017947-IKA001	PETTE Multi	12	10-100	Yes	£663.00	£298.00
0020017948-IKA001	PETTE Multi	12	20-200	Yes	£663.00	£298.00
0020017946-IKA001	PETTE Multi	12	0.5-10	Yes	£663.00	£298.00

DLAB MicroPette Plus and HiPette Manual Pipettes



Code	Model	Max Number of channels	Volume range (min-max)	Autoclavable	List Price	Offer Price
7030301001	MicroPette Plus	1	0.1-2.5	Yes	£74.00	£50.00
7030301004	MicroPette Plus	1	0.5-10	Yes	£74.00	£50.00
7030301005	MicroPette Plus	1	2-20	Yes	£74.00	£50.00
7030301008	MicroPette Plus	1	10-100	Yes	£74.00	£50.00
7030301009	MicroPette Plus	1	20-200	Yes	£74.00	£50.00
7030301014	MicroPette Plus	1	100-1,000	Yes	£74.00	£50.00
7030301016	MicroPette Plus	1	200-1,000	Yes	£74.00	£50.00
7030601001	HiPette	1	0.1-2.5	Yes	£144.00	£86.00
7030601004	HiPette	1	0.5-10	Yes	£144.00	£86.00
7030601005	HiPette	1	2-10	Yes	£144.00	£86.00
7030601008	HiPette	1	10-100	Yes	£144.00	£86.00
7030601009	HiPette	1	20-200	Yes	£144.00	£86.00
7030601014	HiPette	1	100-1,000	Yes	£144.00	£86.00

IN STOCK NOW



Pipetting Controller

DLAB Levo Plus Pipetting Controller



Code	Compatible pipette volumes	Digital display	Speed control	List Price	Offer Price
7061110200	0.1 to 100	Yes	Yes	£230.00	£131.00

IN STOCK NOW



pH and Conductivity Meters



Up to **45 %** Off

Hanna Instruments Edge pH Meter

Reliable, accurate, and easy to use – the Hanna Instruments pH meter is the perfect tool for fast, dependable measurements in any lab or field setting. Compact and intuitive, it delivers precise readings with minimal setup.

Key Features & Benefits:

- Fast, accurate pH measurements
- Simple calibration and operation
- Clear digital display
- Robust, portable design
- Long-lasting electrode performance



Code	Type	Type of probe included	List Price	Offer Price
HI-2020-02	pH	TRIS compatible	£668.00	£365.00

Can measure conductivity and dissolved oxygen with the right probes

IN STOCK NOW



METTLER TOLEDO

Up to **45 %** Off

Mettler SevenDirect pH and Conductivity Meter Kits

30671554 SevenDirect SD20 Kit Includes – pH electrode InLab Expert Pro-ISM and 3×3 calibration buffer sachets pH 4.01 7.00 and 9.21 or 10.01

30671562 SevenDirect SD30 Kit Includes – conductivity sensor InLab 731-ISM and 2×2 conductivity standard sachets

30671567 SevenDirect SD23 Kit Includes – pH electrode InLab expert Pro-ISM conductivity sensor InLab 731-ISM 3×3 pH buffer pH 4.01 7.00 and 9.21 or 10.01 sachets and 2×2 conductivity standards



Code	Type	Type of probe included	List Price	Offer Price
30671554	pH	Standard	£1,452.00	£786.00
30671562	Conductivity	Standard	£1,961.00	£1,088.00
30671567	pH and conductivity	Standard	£2,842.00	£1,624.00

Density Meters



Up to **34%** Off

Mettler Densito Density Meters

Designed to fit comfortably, the lightweight Densito allows for one-handed operations in the lab or on-the-go. Results are temperature-compensated and can be delivered in specific gravity, Brix, API units and others.



Code	Max measuring range (g/cm ³)	Connectivity	List Price	Offer Price
30330857	3	USB	£3,219.00	£2,123.00
30330858	3	USB, RFID and barcode reader	£3,602.00	£2,378.00

Gel Documentation Systems



Up to **29%** Off



Syngene IG5 Gel Documentation System Package

Package includes dedicated Ingenius capture & editing software and GeneTools analysis software (unlimited). Max image size 20 x 24cm, **GE-2020M-E** Transilluminator (20 x 20cm; 302nm) and **GE-CONVERT5** Visible light converter

Code	DNA gels	Protein gels	Chemilluminescence blots	List Price	Offer Price
IG5-Package	Yes	Yes	Yes	£5,817.00	£4,119.00

Bunsen Burners



Up to **29%** Off

Cole-Parmer BB-200-CN Bunsen Burner

Code	Overall height (mm)	Gas type	Controller Included	List Price	Offer Price
36135-02	177	Electric	Yes	£984.00	£695.00



Microscopes



Up to **29%** off

Kern & Sohn OBN Microscopes

The Kern & Sohn OBN microscope series combines professional performance with outstanding flexibility. Designed with ergonomic controls, precision optics and professional Köhler illumination, the OBN range delivers sharp, high-contrast imaging for routine and advanced microscopy work. With modular upgrade options including phase contrast, fluorescence and polarisation capabilities, the OBN series can easily adapt to evolving laboratory requirements.

Key Features & Benefits:

- Infinity plan optical system delivers crisp, high-resolution images
- Professional Köhler illumination ensures excellent light control and image clarity
- Large ergonomic mechanical stage enables precise and comfortable sample positioning
- Expandable modular design allows upgrades for phase contrast, fluorescence and more
- Choice of LED or halogen illumination to suit different laboratory applications



Code	Type	Head type	Max magnification (x)	Camera or tablet included	List Price	Offer Price
OBN135	Standard	Trinocular	1,000	No	£1,398.00	£991.00
OBN132	Standard	Trinocular	1,000	No	£1,432.00	£1,015.00
OBN135C825	Standard	Trinocular	10	Camera	£1,781.00	£1,262.00
OBN132C825	Standard	Trinocular	10	Camera	£1,823.00	£1,293.00
OBN135C832	Standard	Trinocular	10	Camera	£1,909.00	£1,352.00
OBN132C832	Standard	Trinocular	10	Camera	£1,951.00	£1,383.00
OBN132T241	Standard	Trinocular	10	No	£2,335.00	£1,654.00
OBN135T241	Standard	Trinocular	10	Tablet	£2,391.00	£1,624.00
OBN158	Standard	Trinocular	1,000	No	£2,488.00	£1,763.00
OBN159	Standard	Trinocular	100	No	£2,505.00	£1,775.00
OBN135T251	Standard	Trinocular	1,000	Tablet	£2,578.00	£1,751.00
OBN159T251	Standard	Trinocular	1,000	Tablet	£3,688.00	£2,504.00
OBN148	Standard	Trinocular	1,000	No	£5,538.00	£3,923.00
OBN142	Standard	Trinocular	10	No	£5,555.00	£3,935.00
OBN142T251	Standard	Trinocular	1,000	Tablet	£6,887.00	£4,677.00
OBN142C861	Standard	Trinocular	1,000	Camera	£7,865.00	£5,341.00

Kern & Sohn OBT Microscopes

The KERN OBT microscopes combine durable construction, bright LED illumination and precision optics for clear, reliable results. Designed for education and routine laboratory use, the OBT series offers ergonomic operation, smooth focusing and excellent image quality for everyday microscopy applications.



Code	Type	Head type	Max magnification (x)	List Price	Offer Price
OBT101	Standard	Monocular	10	£224.00	£167.00
OBT102	Standard	Monocular	10	£241.00	£179.00
OBT103	Standard	Monocular	10	£271.00	£200.00
OBT104	Standard	Binocular	10	£343.00	£251.00
OBT105	Standard	Monocular	10	£292.00	£215.00
OBT106	Standard	Binocular	10	£360.00	£263.00
OBT231	Standard	LCD display	100	£734.00	£527.00

Optika B-190, B-290 and SLX Higher Education Microscopes

B-190 Series - Advanced student microscopes combining quality, reliability and innovation.

B-290 Series - Durable routine microscopes with brightfield, darkfield, phase contrast and LED fluorescence.

SLX Series - Cordless stereo and stereozoom microscopes for educational, laboratory and industrial applications.



HIGHER EDUCATION MODELS ON SPECIAL OFFER

Code	Type	Head type	Max magnification (x)	List Price	Offer Price
B-190TBPL+	Standard	Binocular	1,000	£1,727.00	£1,119.00
B-290TB	Standard	Binocular	1,000	£1,782.00	£1,138.00
B-293	Standard	Trinocular	1,000	£969.00	£626.00
SLX-1	Stereo	Binocular	40	£436.00	£287.00
SLX-3	Stereo	Trinocular	45	£664.00	£433.00

Optika Higher Education Packages

Explore our Microscope packages designed specifically for higher education.

B-293 Package includes -

Code	List Price	Offer Price
B-293-Package1	£1,291.00	£813.00

Includes **C-B5** camera



B-293

SLX-3 Package includes -

Code	List Price	Offer Price
------	------------	-------------

SLX-3-Package2	£1,340.00	£651.00
-----------------------	------------------	----------------

Includes **M-699** universal adapter for projection lens (trino)

Includes **C-B5** camera



SLX-3

Optika B-380, B-510 and IM-300 Laboratory Hospital Microscopes

B-380 Series - Versatile routine microscopes supporting brightfield, darkfield, phase contrast, fluorescence and polarised light applications, with a durable and stable design for long-term laboratory use.

B-510 Series - Advanced routine microscopes delivering high-performance transmitted light analysis with modern, user-focused design and precision optics.

IM-300 - Ergonomic inverted microscopes designed for comfortable routine observation of living cells and tissue cultures in natural conditions.



LABORATORY HOSPITAL MODELS ON SPECIAL OFFER

Code	Type	Head type	Max magnification (x)	List Price	Offer Price
B-383LD	Standard	Trinocular	1,000	£2,353.00	£1,503.00
B-383PLI	Standard	Trinocular	1,000	£1,262.00	£810.00
B-510BF4K	Standard	Trinocular	1,000	£3,004.00	£1,878.00
IM-300	Inverted	Trinocular	400	£3,061.00	£1,687.00
IM-300D	Inverted	Trinocular	400	£5,402.00	£3,417.00
IM-300LD2	Inverted	Trinocular	400	£8,020.00	£4,996.00

Optika Laboratory Hospital Packages

Explore our Microscope packages designed specifically for hospitals.

B-383PLI Package 1 includes -

Code	List Price	Offer Price
B-383PLI-Package1	£1,694.00	£1,055.00

Includes **M-699** universal adapter for projection lens
Includes **C-B5** camera

B-383PLI Package 2 includes -

Code	List Price	Offer Price
B-383PLI-Package2	£2,298.00	£1,427.00

Includes **M-620.1** 0.5x focusable C-Mount adapter
Includes **C-P6** camera

IM-300LD2 Package includes -

Code	List Price	Offer Price
IM-300LD2-Package	£9,152.00	£5,639.00

Includes **M-620.1** 0.5x focusable C-Mount adapter
Includes **C-P6** camera



B-383PLI



IM-300LD2

Storage Cabinets



Up to **22%** Off

Quality Metal Products Storage Cabinets

Keep hazardous chemicals safely organised and fully compliant with QMP's range of Hazardous Substance, Acid & Alkali and COSHH Security Cupboards. Manufactured from heavy-duty steel with secure locking systems, spill-retaining shelves and liquid-tight sumps, these cupboards provide reliable protection for hazardous materials in laboratories, industrial and commercial environments. Designed for durability and workplace safety, the range helps ensure secure, compliant storage with clear hazard identification and long-lasting performance.



Code	Model	Height (mm)	Width (mm)	Depth (mm)	List Price	Offer Price
CZ703530ZKXX	CoSHH Cupboard	700	350	300	£285.00	£223.00
CA703530ZWXX	Acid and Alkali Cupboard	700	350	300	£306.00	£238.00
CA904646ZWXX	Acid and Alkali Cupboard	900	460	460	£357.00	£277.00
CA189046ZWXX	Acid and Alkali Cupboard	1,800	900	460	£655.00	£503.00
CZ904646ZYXX	Hazardous Cupboard	900	460	460	£325.00	£253.00
CZ909046ZYXX	Hazardous Cupboard	900	900	460	£405.00	£314.00
CZ189046ZYXX	Hazardous Cupboard	1,800	900	460	£582.00	£448.00

Microplate Shakers



Up to **39%** Off

DLab MX-M Microplate Shaker

Code	Model	Max speed (rpm)	List Price	Offer Price
8031202000	MX-M	1,500	£349.00	£213.00



IN STOCK NOW



Rockers

DLAB

Up to **39%**
Off

DLAB SK-R1807-S and SK-D1810-S Rockers

Code	Speed range (rpm)	3D motion	Platform dimensions (w x d)	List Price	Offer Price
8072311203	10-80	No	240 x 215mm	£293.00	£179.00
8032434200	10-80	Yes	240 x 215mm	£205.00	£119.00



IN STOCK NOW



Rotator

DLAB

Up to **39%**
Off

DLAB MX-RL-E Rotator

Code	Max tube diameter (mm)	Max rpm	List Price	Offer Price
8031402000	30	80	£293.00	£179.00



IN STOCK NOW



Roller

DLAB

Up to **39%**
Off

DLAB MX-T6-S+ Roller

Code	Number of positions	Variable speed	Max speed (rpm)	Timer	List Price	Offer Price
8063312000	6	Yes	70	No	£328.00	£200.00



IN STOCK NOW



Aspirators



Up to **39%** Off

DLAB Aspirators

The DLAB SmartVac offers fast, hands-free vacuum filtration with intelligent controls, while the DLAB SafeVac provides safe, adjustable liquid aspiration for everyday lab tasks. Together, they deliver efficient, reliable solutions that simplify workflows and support a cleaner, more user-friendly laboratory environment.



IN STOCK NOW



Code	For hazardous biological waste	For hazardous chemical waste	Flow rate (l/min)	Waste trap volume (l)	List Price	Offer Price
7035300002	Yes	No	1.6	1	£293.00	£179.00
7035100000	Yes	Yes	15	4	£1,094.00	£705.00

Vortexers



Up to **27%** Off

IKA Vortex 2 Vortexer

Code	Vessel size range (ml)	Variable speed	List Price	Offer Price
0025001608-IKA001	0.2-250	Yes	£340.00	£247.00



DLAB MX-S Vortexers



Code	Vessel size range (ml)	Variable speed	Continuous / touch control	List Price	Offer Price
8031102000	0.1-50	Yes	Yes	£205.00	£112.00
8031242100	0.1-50	Yes	Yes	£247.00	£151.00



Up to **44%** Off

IN STOCK NOW



Block Heaters



Up to **25%** Off

Grant Instruments QBD Block Heaters

Grant's QBD Block Heaters deliver consistent, accurate heating with a compact footprint and flexible block options—perfect for routine sample heating and incubation.



Code	Max temp (°C)	Refrigerated	Number of blocks	Blocks included?	List Price	Offer Price
QBD1	130	No	1	No	£840.00	£626.00
QBD2	130	No	2	No	£919.00	£685.00
QBD4	130	No	4	No	£1,119.00	£832.00



Up to **39%** Off

D Lab HB120-S Block Heater

IN STOCK NOW



Code	Max temp (°C)	Refrigerated	Number of blocks	Blocks included?	List Price	Offer Price
5032101200	25	No	1	Yes	£433.00	£264.00

Viscometers



Up to **36%** Off

IKA ROTAVISC Viscometers

Code	Model	Viscosity range (min-max)	Accuracy	Repeatability	List Price	Offer Price
0010013826-IKA001	ROTAVISC lo-vi Complete	1-6,000,000	1	0.2	£3,300.00	£2,102.00
0010013827-IKA001	ROTAVISC me-vi Complete	100-40,000,000	1	0.2	£3,300.00	£2,102.00



Thermal Cyclers



Up to **39%** Off

DLAB TC2000-G1 Thermal Cycler

Reliable, precise, and designed for everyday laboratory use – the DLAB TC2000-G1 Thermal Cycler delivers consistent, accurate temperature control for dependable PCR results. With a compact design and user-friendly interface, it supports efficient workflows in any lab.

Key Features & Benefits:

- Accurate and uniform temperature control for reliable PCR performance
- Intuitive touchscreen interface for easy programming and operation
- Fast heating and cooling rates to reduce run times
- Compact, space-saving design ideal for busy laboratories
- Durable construction for long-term, dependable use



IN STOCK NOW



Code	Type	Capacity 96 well plates	Capacity 0.2ml tubes	List Price	Offer Price
5034102400	Gradient	1	96	£2,776.00	£1,694.00

Cole-Parmer®

Up to **32%** Off

Cole-Parmer PCR-300 Alpha Cyclers

The Cole-Parmer PCR-300 Alpha Cyclers combine flexible multi-block capacity, intuitive touchscreen control and rapid protocol optimisation to deliver reliable, high-throughput PCR performance with minimal training required.



Code	Type	Capacity 384 well plates	Capacity 0.2ml tubes	List Price	Offer Price
93945-02	Gradient	0	96	£4,522.00	£3,066.00
93945-00	Gradient	1	0	£4,522.00	£3,066.00
93945-28	Gradient	2	0	£9,778.00	£6,667.00
93945-20	Gradient	1	96	£9,778.00	£6,667.00
93945-22	Gradient	0	192	£9,824.00	£6,699.00
93945-04	Gradient	3	0	£18,617.00	£12,645.00
93945-06	Gradient	2	96	£18,617.00	£12,645.00
93945-10	Gradient	1	96	£18,705.00	£12,704.00
93945-08	Gradient	4	0	£18,705.00	£12,704.00
93945-12	Gradient	0	96	£18,705.00	£12,704.00

Thermoshakers



Up to **39%** Off

DLab HM100-Pro and HCM100-Pro Thermoshakers

The DLab HM100-Pro and HCM100-Pro Thermoshakers combine precise temperature control with powerful mixing performance. Delivering fast, reliable sample preparation for molecular biology, enzymatic reactions and PCR applications, while programmable settings, automatic block recognition and compatibility with a wide range of tubes maximise laboratory efficiency and flexibility.



Code	Model	Temperature range (°C)	Shaking speed range (rpm)	Maximum orbit (mm)	List Price	Offer Price
5063104300	HM100-Pro	0-100	200-1,500	3	£849.00	£518.00
5062143300	HCM100-Pro	26-100	200-1,500	3	£1,080.00	£659.00

IN STOCK NOW



Up to **28%** Off

IKA MATRIX Orbital Delta Thermoshaker

The IKA MATRIX Orbital Delta thermoshaker delivers reliable mixing and precise temperature control in one compact unit, ensuring fast, contamination-free sample preparation for DNA/RNA, ELISA and pharmaceutical applications with even the smallest sample volumes.

Key Features & Benefits:

- **High-speed orbital shaking up to 3000 rpm** – provides rapid, consistent sample mixing for improved reproducibility.
- **Precise heating from 25°C to 100°C** – ideal for sensitive biological and molecular workflows requiring accurate temperature stability.
- **User-friendly display and intuitive controls** – simplifies operation and reduces setup time in busy laboratories.
- **Compact, robust aluminium housing** – ensures excellent stability, durability and reliable everyday performance.



Code	Temperature range (°C)	Shaking speed range (rpm)	Maximum orbit (mm)	List Price	Offer Price
0020100517-IKA001	1-100	0-3,000	3	£2,565.00	£1,856.00

Homogenisers



Up to **39%** Off

DLab D-160 Package 2 Homogeniser

Code	Motor (watts)	Volume range (min-max)	List Price	Offer Price
8034121000	160	1-250	£1,185.00	£723.00

Package includes **D160** drive, **H 370** stand and **DS-160/10** shaft

IN STOCK NOW



Cryostats



Up to **9%** Off

Bright Instruments OTF7000/LS Cryostat

The OTF7000/LS cryostat delivers highly stable temperature control and precision sectioning, combined with UV disinfection, intuitive touchscreen operation and efficient specimen cooling for reliable, reproducible results in demanding clinical and research applications.

Code	Touch Screen	Maximum specimen size (mm)	List Price	Offer Price
OTF7000/LS-001	Yes	70	£34,028.00	£30,965.00



Can't Find What You're Looking For?

Our brochure highlights just a selection of the products and offers available from WolfLabs. If you can't find the equipment you need, require further information, or would like a tailored quotation, our friendly team is here to help.

Visit our Contact Us page: www.wolf labs.co.uk/contact-us

sales@wolf labs.co.uk

01759 301142

Get in touch today and let us help you find the right solution for your laboratory.